



**PARADIGM™**

**Detection Platform**

*Multimode Plate Detection*

葉佳振

*Beckman Coulter, Taiwan*

# Paradigm

- Absorbance
- Fluorescence Intensity
- Time Resolved Fluorescence
- Fluorescence Resonance Energy Transfer (FRET)
- Fluorescence Polarization
- HTRF
- Luminescence

# Paradigm



- Base Instrument

- Detection Cartridge

- Top Read (6 positions)
- Bottom Read (6 positions)



# FRONT:

## Base Instrument Interface



Touch buttons for plate and detection cartridge access

Status light

Full front door access for flexible positioning for robotic integration

Plate height sensor

4 plate orientations





# Base Instrument




Lead the Way  
Together in AP



## HIGHLIGHTS

- User setup and install
- Z-Height Adjustment (top & bottom)
- Plate Height Sensing
- Intelligent Cartridge Recognition
- Dual PMT System
- CO<sub>2</sub> Atmosphere Option
- Shaking – Linear & orbital
- Temperature Control
  - Ambient to +45
- 6 ~ 1,536 well Plates
- USB or com port connection

- 獨立式卡匣光學偵測模組





# Remove Cartridge





# Install Cartridge



- **Existing Detection Cartridges**
  - **Absorbance (ABS)**
  - **Multimode (MULTI)**
    - **FI, LUM, TRF**
  - **Fluorescence Intensity (FI)**
  - **Fluorescence Polarization (FP)**
  - **Time Resolved Fluorescence (TRF)**
  - **Luminescence (LUM)**
  - **HTRF® - Cisbio**





# Absorbance (ABS) Detection Cartridge



## PARADIGM™ Detection Platform Absorbance Detection Cartridge

### DESCRIPTION

Part Number.....A41575  
Mounting.....Top Read (Photodiode Detector) and  
Bottom Read (Monochromatic Illumination)  
Cartridge Size.....Top - 1 Position  
Bottom - 2 Position

### TYPICAL PERFORMANCE

Wavelength Range.....230 - 1000 nm  
Wavelength Selection.....Monochromator (1.0 nm increments)  
Wavelength Bandwidth.....4.0 nm  
Wavelength Accuracy..... $\pm 1$  nm  
Wavelength Repeatability..... $\pm 0.5$  nm  
Photometric Range.....0 - 3.5 OD  
Photometric Resolution.....0.0001 OD  
Photometric Accuracy..... $< \pm 0.01$  OD  $\pm 2\%$ , 0-20D @ 405 nm  
Photometric Precision..... $< \pm 0.01$  OD  $\pm 0.5\%$ , 0-20D @ 405 nm



Covering the UV/Visible wavelength range necessary to perform a variety of direct measurements of protein and nucleic acid applications fast and easy, the state of the art monochromatic illumination design quickly adjusts to address a variety of ELISA and immunoassay wavelengths.



# Multimode (MULTI) Detection Cartridge



Lead the Way  
Together in AP



PARADIGM™ Detection Platform  
Multimode Detection Cartridge

## DESCRIPTION

Part Number.....	A41576
Mounting.....	Top or Bottom Read
Cartridge Size.....	3 Position
Light Source.....	High-Power LED
Wavelength-Tuned Excitation/Emission Range .....	360/35 nm - 465/35 nm 485/20 nm - 535/25 nm 535/35 nm - 595/35 nm 585/10 nm - 635/35 nm 370/80 nm - 616/10 nm

Optimized for 96- and 384-density assays, the Multimode cartridge provides a broad spectrum of single emission read capabilities, including:

1. Fluorescence Intensity (FI)
2. Time Resolved Fluorescence (TRF)
3. Glow Luminescence (LUM)

## TYPICAL APPLICATIONS\*\*

Proliferation/Viability  
Nucleic Acid Quantitation  
cAMP Quantitation  
GPCR  
Immunoassay  
Ion Channel  
Reporter

## TYPICAL FLUOROPHORES\*\*

Coumarin  
Fluorescein  
Texas Red  
Rhodamine  
Europium (TRF only)



# Fluorescence Intensity (FI) Detection Cartridge



Lead the Way  
Together in AP



Optimized for ultra fast and sensitive detection in 384 and 1536 formats, it's dual emission design performs both single and dual emission assays in a single read of the plate. Making *single emission* **Fluorescence Intensity (FI)** and *dual emission* **Fluorescence Resonance Energy Transfer (FRET)** measurements simple and straightforward.

## PARADIGM™ Detection Platform Fluorescence Intensity Detection Cartridge

### DESCRIPTION

Part Number.....	A41577 & A41578
Mounting.....	Top Read or Bottom Read
Cartridge Size.....	1 Position
Light Source.....	High Powered LED
A41577:	
Wavelength Tuned Excitation Range.....	360/35 nm
Wavelength Tuned Emission Range.....	465/35 nm
	535/25 nm
A41578:	
Wavelength Tuned Excitation Range.....	485/20 nm
Wavelength Tuned Emission Range.....	535/25 nm
	595/35 nm

TYPICAL APPLICATIONS**	TYPICAL FLUOROPHORES**	TYPICAL FLUOROPHORES**
Both A41577 & A41578: Nucleic Acid Quantitation Molecular Interaction A41578 only: Ion Channel Protein Quantitation Reporter Apoptosis Proliferation/Viability	A41577: Coumarin Hoechst DAPI AMCA ANS Coumarin-fluorescein (FRET)	A41578: Fluorescein Oligreen Picogreen Ribogreen Fluo-3,4 GFP Fluorescein-Rhodamine (FRET) GFP-RFP (FRET)



# Fluorescence Polarization (FP) Detection Cartridge



Lead the Way  
Together in AP

Optimized for ultra fast and sensitive detection in 384 and 1536 formats, it's dual emission optical design allows the PARADIGM to simultaneously collect all necessary parallel and perpendicular emission data necessary for **fluorescence polarization** detection.

## PARADIGM™ Detection Platform Fluorescence Polarization Detection Cartridge

### DESCRIPTION

Part Number.....	A41581 or A41582
Mounting.....	Top Read or Bottom Read
Cartridge Size.....	1 Position
Light Source.....	High-Power LED
A41581	
Wavelength-Tuned Excitation Range .....	485/20 nm
Wavelength-Tuned Emission Range.....	//535/25 nm
	⊥ 535/25 nm
A41582	
Wavelength-Tuned Excitation Range .....	535/25 nm
Wavelength-Tuned Emission Range.....	//595/35 nm
	⊥ 595/35 nm

### TYPICAL APPLICATIONS\*\*

**A41581 & A41582:**  
cAMP Quantitation  
Genotyping (AcycloPrime)  
Immunoassay  
Molecular Interaction  
Nuclear Receptor

### TYPICAL FLUOROPHORES\*\*

**A41581:**  
Fluorescein  
FP - green

**A41582:**  
Rhodamine  
FP - red



# Time Resolved Fluorescence (TRF) Detection Cartridge



Optimized for ultra fast and sensitive detection in 384 and 1536 formats, it's single or dual emission design performs both single and dual emission assays in a single read of the plate. Making *single emission* or *dual emission* **time resolved fluorescence (TRF)** simple and fast.

## PARADIGM™ Detection Platform Time Resolved Fluorescence Detection Cartridge

### DESCRIPTION

Part Number.....	A41579
Mounting.....	Top Read or Bottom Read
Cartridge Size.....	2 Position
Light Source.....	High Powered LED
Wavelength Tuned Excitation.....	370/80 nm
Wavelength Tuned Emission.....	616/10 nm 642/10 nm

### TYPICAL APPLICATIONS\*\*

cAMP Quantitation  
Immunoassays  
GTP Binding  
Apoptosis  
GPCR Ligand Binding

### TYPICAL FLUOROPHORES\*\*

Europium  
Samarium



# Typical Fluorophores




Lead the Way  
Together in AP


Coumarin	Fluorescein	Rhodamine	Europium
Hoechst	OliGreen	Texas Red	Europium Chelates
DAPI	PicoGreen	TAMRA	Samarium Chelates
AMCA	RiboGreen	AlexaFluor®532-555	Coumarin-Fluorescein(FRET)
ANS	Fluo -- 3, 4	Cy3	Fluorescein-Rhodamine(FRET)
AlexaFluor®350	GFP	R-Phycoerythrin	GFP-RFP(FRET)
4-MU	CellTiter Blue	Resorufin	
	AlamarBlue	AlexaFluor®594	
	EFC	Spectrum Red	
	FITC		
	BODIPY FL		
	Calcein		
	Cy2		
	OregonGreen		
	Sytox Green		

# Optical system - Slit

- 單光柵可將光的純度提高1000倍
- 雙光柵可將光的純度提高 $10^6$ 倍, 但相對的會造成光能量降低, 影響靈敏度



# Monochromator




AGif - UNREGISTERED



Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



### Fluorophore

- |         |  |  |
|---------|--|--|
| 1: DAPI | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

### Excitation

Laser (nm): None

Filter / Bandpass

360 / 5

### Emission

Filter / Bandpass

465 / 12

Toggle Excitation Display

Toggle Emission Display

Clear All


### Legend

- Absorption and fluorescence emission spectra of DAPI bound to DNA.

Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



## Fluorophore

- |         |  |  |
|---------|--|--|
| 1: DAPI | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

## Excitation

Laser (nm): None

Filter / Bandpass

360 / 35

## Emission

Filter / Bandpass

465 / 35

Toggle Excitation Display

Toggle Emission Display

Clear All


## Legend

- Absorption and fluorescence emission spectra of DAPI bound to DNA.

Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



### Legend

Absorption and fluorescence emission spectra of Hoechst 33258 bound to DNA.

### Fluorophore

- |                  |  |  |
|------------------|--|--|
| 1: Hoechst 33258 | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

### Excitation

Laser (nm): None

Filter / Bandpass

360 / 5

### Emission

Filter / Bandpass

465 / 12

Toggle Excitation Display


Toggle Emission Display

Clear All

Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



### Legend

Absorption and fluorescence emission spectra of Hoechst 33258 bound to DNA.

### Fluorophore

- |                  |  |  |
|------------------|--|--|
| 1: Hoechst 33258 | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None          | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

### Excitation

Laser (nm): None

Filter / Bandpass

360 / 35

### Emission

Filter / Bandpass

465 / 35

Toggle Excitation Display


Toggle Emission Display

Clear All

Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



### Legend

Absorption and fluorescence emission spectra of Rhodamine Green™ carboxylic acid (5-(and-6)-carboxyrhodamine 110) in pH 7.0 buffer.

### Fluorophore

- |                    |  |  |
|--------------------|--|--|
| 1: Rhodamine Green | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

### Excitation

Laser (nm): None

Filter / Bandpass

485 / 5

### Emission

Filter / Bandpass

535 / 12

Toggle Excitation Display


Toggle Emission Display

Clear All

Now you can plot and compare spectra and check the spectral compatibility for many fluorophores offered by Molecular Probes.

The Spectra Viewer can most easily be printed by capturing a screen-shot and printing the resulting image file.

For printing instructions or to answer other questions you have, see our [User Guide](#).



### Legend

Absorption and fluorescence emission spectra of Rhodamine Green™ carboxylic acid (5-(and-6)-carboxyrhodamine 110) in pH 7.0 buffer.

### Fluorophore

- |                    |  |  |
|--------------------|--|--|
| 1: Rhodamine Green | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 2: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 3: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 4: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |
| 5: None            | <input checked="" type="checkbox"/> ex | <input checked="" type="checkbox"/> em |

### Excitation

Laser (nm): None

Filter / Bandpass

485 / 20

### Emission

Filter / Bandpass

535 / 25


Toggle Excitation Display

Toggle Emission Display


Clear All

## Excitation Maximum


Excitation ranges and maxima



## Emission Range



## Fluorescence Emission



# Paradigm vs T

	Paradigm	Va--
Wavelength Selection	Filter	Monochromator
Bandwidth	20~35nm	Ex 5nm Em12nm
Sensitivity	5 amol/well	120 amol/well



# Luminescence (LUM) *Detection Cartridge*



Utilizing a unique, well isolating design that restricts the level of cross talk from surrounding wells. The resulting performance is ideal for chemiluminescence glow applications.

## PARADIGM™ Detection Platform Luminescence Detection Cartridge

### DESCRIPTION

Part Number ..... A41584  
Mounting ..... Top Read or Bottom Read  
Cartridge Size ..... 1 Position  
Emission Range ..... Visible to 800 nm

### TYPICAL APPLICATIONS\*\*

Apoptosis  
cAMP Quantitation  
GPCR Ligand Binding  
GTP Binding  
Immunoassay



# 100% User Configurable



- Configure & Expand a system's capabilities in < 2minutes
- Detection Cartridge Positions
  - 6 – Top
  - 6 – Bottom

“PLUG and DETECT”

Get what you need, when you need it . . .





- Easy integrated with Automation
- Expand throughput

*PARADIGM Integrated with a Biomek NX<sup>P</sup> S8*


# What is On the Fly?

## Stop & Go



*12 discrete movements required to measure  
the row of a 96-well plate*

## On the Fly



*1 discrete movements is performed while 12  
discrete measurements are taken*



# Faster Detection!



Added the “On The Fly” Detection to

- Optimized for Speed
- Optimized for Performance

	<u>Stop &amp; Go</u>		<u>Speed</u>		<u>Performance</u>	
	Typical Performance (fmol/well)	Throughput (seconds)	Typical Performance (fmol/well)	Throughput (seconds)	Typical Performance (fmol/well)	Throughput (seconds)
96-well	0.15	35	--	20	0.3	30
384-well	0.1	120	--	30	0.2	50
1536-well	0.03	240	--	55	0.1	90

*NOTE: Data is representative of fluorescence intensity detection cartridge, A41577, with fluorescein*

# Paradigm

- **Multimode Analysis Software**
  1. Drag & Drop protocol work flow development
    - Multiplexing of read modes in a single protocol
  2. Auto update function ( protocols,labware, application notes )
  3. 3D Data display
  4. Labware Optimization
  5. Read Height Optimization

# Drag – Drop Workflow



Create Protocol Training Protocol



General Settings  
Technique Type  
Labware Selection  
Layout Settings  
**Method Selection** >  
Data Reduction Page  
Transformation  
Output Settings

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

Single Kinetic Area Scan Wavelength Scan

Select Method

- [-] Preparation
  - [-] Wait
  - [-] Set Temperature
  - [-] Shake
  - [-] Eject
  - [-] Load
  - [-] Pause
- [+] AbsMono
- [-] FP-FLUO
  - FP - fluorescein(fly)
  - FP - fluorescein(400ms)
- [+] FP-RHOD
  - FI - fluorescein+rhodamine(400ms)
  - FI - fluorescein+rhodamine(fly)
- [+] FI-COFL
- [+] HTRF
- [+] MULTI

☒ Cartridge View  
☐ Show only runnable Methods



Estimated Time 00:08:45

Group1

- [-] Temperature
- FP - fluorescein(400ms)
- Shake
- Area Scan
- FP - fluorescein(fly)

Report Options

### Temperature Info

Set Temperature 37  
Wait for Temperature False

### Set Temperature

Enter the desired microplate chamber temperature in degrees Celsius. (0 turns it off).

\$229

Cancel < Back Next >



Create Protocol NewProtocol 1

General Settings  
Technique Type  
Labware Selection  
Layout Settings  
**Method Selection**  
Data Reduction Page  
Output Settings

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

Single **Kinetic** Area Scan Wavelength Scan

Select Method

Preparation

- Wait
- Set Temperature
- Shake
- Eject
- Load
- Pause
- ABS-MONO
  - Abs\_595nm
  - Abs\_mono\_580nm
  - DNA 2 @977
  - DNA 2 @900
  - DNA 2 @320
  - DNA 2 @280
  - DNA 2 @260
  - Abs method
  - x\_Abs\_Mono\_405nm
  - x\_Abs\_Mono\_280nm
  - x\_Abs\_Mono\_260nm

- ☒ Cartridge View  
☒ Show only runnable Methods



Estimated Time 00:16:51

Group1

- Kinetic
- Abs\_595nm

Report Options

**Kinetic Info**

Kinetic Cycles **50**  
Kinetic Interval **20**

**Kinetic Cycles**

Choose the number of measurement cycles to perform.

Cancel < Back Next >

# Area scan



Create Protocol NewProtocol 1

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

General Settings  
Technique Type  
Labware Selection  
Layout Settings  
**Method Selection**  
Data Reduction Page  
Output Settings

Single Kinetic **Area Scan** Wavelength Scan

Select Method

- Preparation
  - Wait
  - Set Temperature
  - Shake
  - Eject
  - Load
  - Pause
- ABS-MONO
  - Abs\_595nm
  - Abs\_mono\_580nm
  - DNA 2 @977
  - DNA 2 @900
  - DNA 2 @320
  - DNA 2 @280
  - DNA 2 @260
  - Abs method
  - x\_Abs\_Mono\_405nm
  - x\_Abs\_Mono\_280nm
  - x\_Abs\_Mono\_260nm

Estimated Time 00:00:36

Group1

- Area Scan
  - Abs\_595nm

Report Options

Area Scan Info

Area Scan Points 2x2[0.225mm]

**Area Scan Points**  
The number of scan points in each well and scan resolution in mm.

Cancel < Back Next >

Windows taskbar: 開始, 新資..., 10 OS..., Micros..., Multim..., 未命名..., CE, 96%, 上午 10:25

# Area scan



Create Protocol NewProtocol 1

**Method Selection**

Available detection and preparation methods are listed on the left. Click the Add button to add a method to the protocol, or select the method to the protocol, or select the method to the protocol.

**Area Scan Selection Editor**

Select the scan resolution and measurement points to read. The type of linear scan that may be configured is determined by which reading direction is configured in Layout Settings. When Read by row is selected, a horizontal linear scan may be configured. When Read by column or read by well is selected, a vertical linear scan may be configured.

Resolution (mm): 0.95625

Selected Scan Points: Points X: 6, Points Y: 6

Cartridge View: ☒ Show only runnable Methods: ☒

Area Scan Info: 2x2[0.225mm]

Scan Points: Number of scan points in each well and scan resolution in mm.

OK Cancel

Cancel < Back Next >

Windows taskbar: 開始 新資... 10 OS... Micros... Multim... 未命名... CE 96% 上午 10:27

# Area scan



Create Protocol NewProtocol 1

**Method Selection**

Available detection and preparation methods are listed below. Select a method and click the Add button.

Single Kinetic Area Scan Well Scan

Select Method

- Preparation
  - Wait
  - Set Temperature
  - Shake
  - Eject
  - Load
  - Pause
- ABS-MONO
  - Abs\_595nm
  - Abs\_mono\_580nm
  - DNA 2 @977
  - DNA 2 @900
  - DNA 2 @320
  - DNA 2 @280
  - DNA 2 @260
  - Abs method
  - x\_Abs\_Mono\_405
  - x\_Abs\_Mono\_280
  - x\_Abs\_Mono\_260

☒ Cartridge View  
☒ Show only runnable Methods

**Area Scan Selection Editor**

Select the scan resolution and measurement points to read. The type of linear scan that may be configured is determined by which reading direction is configured in Layout Settings. When Read by row is selected, a horizontal linear scan may be configured. When Read by column or read by well is selected, a vertical linear scan may be configured.

Resolution (mm): 0.22500

Selected Scan Points: Points X: 20, Points Y: 20

OK Cancel

**Area Scan Info**

Area Scan Points: 2x2[0.225mm]

**Report Options**

**Scan Points**

Number of scan points in each well and scan resolution in mm.

Cancel < Back Next >

Windows taskbar: 開始, 新資..., 10 OS..., Micros..., Multim..., 未命名..., CE, 96%, 上午 10:26

# Wavelength Scan



Create Protocol NewProtocol 1

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

General Settings  
Technique Type  
Labware Selection  
Layout Settings  
**Method Selection**  
Data Reduction Page  
Output Settings

Single Kinetic Area Scan **Wavelength Scan**

Select Method

- Preparation
  - Wait
  - Set Temperature
  - Shake
  - Eject
  - Load
  - Pause
- ABS-MONO
  - scan
  - x\_Abs\_Wavelength\_Scan**

Estimated Time 00:00:43

Group1

- Wavelength Scan
  - x\_Abs\_Wavelength\_Scan**

Report Options

### Method Info

Minimum Wavelength	230
Maximum Wavelength	1000
Number of Ranges	1
Wavelength Increme	1
Start Wavelength	230
End Wavelength	1000

**Number of Ranges**  
Number of ranges that should be measured.

☒ Cartridge View  
☒ Show only runnable Methods

Cancel < Back Next >

Windows taskbar: 開始, 新資..., 10 OS ..., Micros..., Multim..., 未命名..., CE, 96%, 上午 10:28

# Wavelength Scan



Create Protocol NewProtocol 1

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

General Settings  
Technique Type  
Labware Selection  
Layout Settings  
**Method Selection**  
Data Reduction Page  
Output Settings

Single Kinetic Area Scan **Wavelength Scan**

Select Method

- Preparation
  - Wait
  - Set Temperature
  - Shake
  - Eject
  - Load
  - Pause
- ABS-MONO
  - scan
  - x\_Abs\_Wavelength\_Scan

Estimated Time 00:00:38

Group1

- Wavelength Scan
  - x\_Abs\_Wavelength\_Scan

Report Options

**Method Info**

Minimum Wavelength	230
Maximum Wavelength	1000
Number of Ranges	3
Wavelength Increment	2
Start Wavelength	230
End Wavelength	380
Range 2 Start Wave	490
Range 2 End Wave	600
Range 3 Start Wave	750
Range 3 End Wave	1000

**Number of Ranges**  
Number of ranges that should be measured.

Cartridge View ☒  
Show only runnable Methods ☒

Cancel < Back Next >

Windows Taskbar: 開始, 新資..., 10 OS..., Micros..., Multim..., 未命名..., CE, 96%, 上午 10:30

# Grouping Multiple Reads

Create Protocol Training Protocol

Way in AP

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

Report Options

Single Kinetic Area Scan Wavelength Scan

Select Method

- Preparation
  - Wait
  - Set Temperature
  - Shake
  - Eject
  - Load
  - Pause
- AbsMono
  - NewMethod 3
  - DNA 2 @900
  - DNA 2 @280
  - DNA 2 @260
  - NewMethod 2
  - NewMethod 1
  - DNA@320
- FP-FLUO
- FP-RHOD
- FI-FLRH
- FI-COFL
- HTRF
- MULTI

☒ Cartridge View

☐ Show only runnable Methods



Estimated Time	00:07:50	Method Info
Group0		
	DNA 2 @260	Excitation
Group1		
	DNA 2 @260	

Pre-read with Buffer

Read with Buffer & Sample

*Enables the expanded integration capabilities as well as improved multi-step reaction analysis*

Method Info

S229

Cancel

< Back

Next >

# Data Reduction



Create Protocol Training Protocol

## Data Reduction Page

Press F1 for more information about data reduction functions and formulas.

Report Options

[Add new Pass](#)

Group0

A = DNA 2 @260

Group1

B = DNA 2 @260

First Pass

REDUCTION\_A1

Formula

Formula

B-A

Name of Data Sample - PreRead

Name of Units OD

Notes

This represents a "pre-read" of the buffer, sample addition and then read of the "sample" plate. Reduction\_A1 is the subtraction of sample - buffer

Apply Formula for Wells with Category

- ☐ Blank
- ☐ Control
- ☐ Empty
- ☐ Negative Control
- ☐ Positive Control
- ☒ Sample
- ☐ Standard

[Add new Item](#)

S229


Cancel

< Back


Next >

Excel Formula Syntax

# 3D Display





# 3D Display






# Auto Update

(Protocols, Labware, Application Notes)





# Labware Optimization



# Read Height Optimization (FI,FP,TRF,LUM)

Read Height Optimization Wizard

Prepare the Labware  
Optimize  
**Optimization Complete**

Optimization is complete. Click Save to save the optimized read height.

Optimized Read Height: 5.36 mm      Custom Read Height: 5.35 mm

Cancel   Back   Save

Method: Luminescence

# Run Protocol

Run Protocol Invitrogen Fluorescence Assays

Assay Selection Page  
Prepare to Run >  
Run Protocol

## Prepare to Run

Click Run to run the protocol.

Result Name: 20060614-141235

# Samples to measure: 1

Eject Plate Carrier

Close Plate Carrier

Run the selected protocol

Select Plate Orientation (click on Plate to view Layout)

☒ Landscape (Optimized)  
☐ Portrait  
☐ Opposite Landscape  
☐ Opposite Portrait

Instrument

Optimization Features

☒ Do Labware Optimization before Run. Labware Optimization should be performed every Month to get best performance of your instrument. Last optimization was performed on 14.06.2006.

☒ Do Read Height Adjust before Run. Read Height adjustment should be performed on each volume change to get best performance of your instrument. Currently 1,0 mm above Plate performed on 14.06.2006.

Cancel < Back Run >

Input # of Samples

Plate Layout

# Paradigm

- Create method and protocol

Multimode Analysis Software

File Actions Help

Create

Edit

Copy

Delete

Protocols

Detection Methods

Results

Labware

Instruments

Trash

Detection Method Selection List

Name	Measurement Technique	Cartridge	Created	Last Edited	Enabled
Lum_1000	Luminescence	MULTI (s/n 1010)	7/25/2007 2:49:40 PM	7/25/2007 2:50:44 PM	
x_Abs_Mono_382nm	Absorbance	ABS-MONO (s/n 1009)	7/25/2007 1:50:23 PM	7/25/2007 1:53:15 PM	
x_Abs_Wavelength_Scan	Absorbance	ABS-MONO (s/n 1009)	5/16/2007 7:26:46 PM	5/16/2007 7:27:20 PM	
x_Abs_Mono_405nm	Absorbance	ABS-MONO (s/n 1009)	5/16/2007 7:25:28 PM	5/16/2007 7:26:39 PM	
x_Abs_Mono_280nm	Absorbance	ABS-MONO (s/n 1009)	5/16/2007 7:25:28 PM	5/16/2007 7:26:24 PM	
x_Abs_Mono_260nm	Absorbance	ABS-MONO (s/n 1009)	5/16/2007 7:25:28 PM	5/16/2007 7:26:06 PM	
x_Multi_Texas Red Int Top 140m	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 7:56:56 PM	5/11/2007 11:40:58 PM	
x_Multi_Rhodamine Int Top 140	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 6:03:25 PM	5/11/2007 11:40:43 PM	
x_Multi_Luminescence 140ms	Luminescence	MULTI (s/n 1010)	5/11/2007 4:59:02 PM	5/11/2007 11:26:22 PM	
x_Multi_Europium Chelate 125pl	Time Resolved Fluorescence	MULTI (s/n 1010)	5/11/2007 5:13:33 PM	5/11/2007 9:33:37 PM	
x_Multi_Coumarin Int Top 140m	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 6:07:26 PM	5/11/2007 9:31:14 PM	
x_Multi_Fluorescein Int Top 140	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 5:21:27 PM	5/11/2007 9:30:53 PM	
x_Multi_Coumarin Int Top 400m	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 6:07:26 PM	5/11/2007 8:02:03 PM	
x_Multi_Texas Red Int Top 400m	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 7:56:56 PM	5/11/2007 7:57:59 PM	
x_Multi_Rhodamine Int Top 400	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 6:03:25 PM	5/11/2007 6:04:48 PM	
x_FI_RET Coum-Fluor 140ms	FRET	FI-COFL (s/n 1009)	5/11/2007 5:57:47 PM	5/11/2007 5:58:23 PM	
x_FI_Coumarin Int Top 250ms	Fluorescence Intensity Top	FI-COFL (s/n 1009)	5/11/2007 5:54:06 PM	5/11/2007 5:55:33 PM	
x_FI_Coumarin Int Top 140ms	Fluorescence Intensity Top	FI-COFL (s/n 1009)	5/11/2007 5:54:06 PM	5/11/2007 5:55:21 PM	
x_FI_Coumarin Int Top fly perf	Fluorescence Intensity Top	FI-COFL (s/n 1009)	5/11/2007 5:54:06 PM	5/11/2007 5:54:24 PM	
x_Multi_Fluorescein Int Top 400	Fluorescence Intensity Top	MULTI (s/n 1010)	5/11/2007 5:21:27 PM	5/11/2007 5:22:27 PM	
x_Multi_Europium Chelate 1000p	Time Resolved Fluorescence	MULTI (s/n 1010)	5/11/2007 5:13:33 PM	5/11/2007 5:16:00 PM	
x_Multi_Europium Chelate 400pl	Time Resolved Fluorescence	MULTI (s/n 1010)	5/11/2007 5:13:33 PM	5/11/2007 5:15:39 PM	
x_Multi_Luminescence 1000ms	Luminescence	MULTI (s/n 1010)	5/11/2007 4:59:02 PM	5/11/2007 5:12:57 PM	

Parameters

Excitation Wavelength (nm)485

Emission Wavelength (nm)535

Integration Time400.00 ms

Hide Preview

LTTW03A0108

PARADIGM 5/N:1012

Top


Bottom


start


Multimode Analysis S...


untitled - Paint

3:21 PM









# Protocol

- Protocol包含:
  - 偵測模式( Abs, FI, FP, TRF, Lum, HTRF )
  - 是否需要控溫, Shake, Pause
  - 待測微孔盤的形式( 96, 384, 1536 )
  - 樣品Layout( Control, Standard, Blank, Positive control, Negative control, Sample, Dilute... )
  - 資料計算, 輸出列印形式

Multimode Analysis Software

File Actions Help

Create Run Edit Copy Print Delete

Protocols

Detection Methods

Results

Labware

Instruments

Trash

## Protocol Selection List

Name	Application Type	Created	Last Edited	Enabled
Area Scan	Analysis	7/23/2007 5:25:25 PM	7/25/2007 6:26:33 PM	
DNA Ratio	Analysis	7/25/2007 5:48:30 PM	7/25/2007 6:04:05 PM	
Lum_1000	Analysis	7/25/2007 2:50:53 PM	7/25/2007 3:35:07 PM	
SCAN 350 TO 1000	Analysis	7/23/2007 4:43:46 PM	7/25/2007 1:35:11 PM	
x_Abs_405nm_384well	Analysis	5/14/2007 1:27:25 PM	5/16/2007 7:30:52 PM	
x_Abs_405nm_96well	Analysis	5/14/2007 1:27:25 PM	7/25/2007 4:53:46 PM	
x_Abs_Gen_96well	Analysis	5/14/2007 1:27:25 PM	5/16/2007 7:29:39 PM	
x_Abs_Gen_96well_w_scan	Analysis	5/14/2007 1:27:25 PM	5/16/2007 7:28:54 PM	
x_Abs_Mono_382nm_96Well	Analysis	7/25/2007 1:53:59 PM	7/27/2007 12:07:41 PM	
x_FI_Coumarin Int Top 1536w 1.5min	Analysis	5/8/2007 5:02:32 PM	5/12/2007 10:59:52 PM	
x_FI_Coumarin Int Top 384w 1min	Analysis	5/8/2007 5:02:32 PM	5/12/2007 10:59:29 PM	
x_FI_Coumarin Int Top 384w 2min	Analysis	5/8/2007 5:02:32 PM	5/12/2007 10:59:09 PM	
x_FI_Coumarin Int Top 96w 1min	Analysis	5/8/2007 5:02:32 PM	5/12/2007 10:58:44 PM	
x_FI_FRET_Coum-Fluor Int Top 96w 1min	Analysis	5/8/2007 5:02:32 PM	5/21/2007 3:03:15 PM	
x_Multi_Coumarin Int Top 384well 2min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:29:56 PM	
x_Multi_Coumarin Int Top 96well 1.5min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:29:08 PM	
x_Multi_Europium Chelate 384well 2min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:28:47 PM	
x_Multi_Europium Chelate 384well 4min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:28:25 PM	
x_Multi_Europium Chelate 96well 2.5min	Analysis	5/11/2007 8:36:55 PM	7/27/2007 5:36:13 PM	
x_Multi_Fluorescein Int Top 384well 2min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:26:45 PM	
x_Multi_Fluorescein Int Top 384well 4min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:26:24 PM	
x_Multi_Fluorescein Int Top 96well 1.5min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:26:08 PM	
x_Multi_Luminescence 384well 2min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:25:06 PM	
x_Multi_Luminescence 384well 4min	Analysis	5/11/2007 8:36:55 PM	5/12/2007 10:24:46 PM	

**Parameters**

Application Analysis

Labware x\_FI\_Coum\_Greiner 96 black std\_

State Normal

Instrument Unknown

Notes quick start protocol; make sure the labware definition matches your plate height; add standard to corner wells to perform labware and

Hide Preview

LTTW03A0108 PARADIGM 5/N:1012

Top Bottom

start

Multimode Analysis S...

Paradigm Instrument ...

3:20 PM

## General Settings

Please enter a name and notes for this protocol.

Protocol name

Date Created Thursday, August 16, 2007

Date Edited Thursday, August 16, 2007

Date last run Thursday, August 16, 2007

Notes

### Analysis Options

- ☐ Variables
- ☐ Transformation
- ☒ Concentration
- ☒ Cutoff
- ☒ Validation

Cancel

Next >

General Settings

Technique Type &gt;

Labware Selection

Layout Settings

Method Selection

Data Reduction Page

Concentration

Cutoff







Validation

Output Settings

## Technique Type

Select the desired technique type from the list below.

Technique Type

- ☒  Absorbance
- ☐  Luminescence
- ☐  FRET
- ☐  Fluorescence Intensity Top
- ☐  Fluorescence Intensity Bottom
- ☐  Time Resolved Fluorescence

Cancel

&lt; Back

Next &gt;

## Labware Selection

Select the desired labware type from the list below.

Type of Labware

Name	Microplate Format
Standard 384	384
Standard 1536	1536
x_Abs_Greiner 384 VIS clear std	384
x_Abs_Greiner 96 UV clear std	96
Standard 96	96
x_Abs_Greiner 96 VIS clear std	96
ChemLib 1536	1536

Cancel



Back

Next



## Layout Settings

Import Layout

Report Options

Type: Sample

Index: 73

Filling

☐ Vertical☒ Horizontal

Flow

☐ Constant☒ Incremental

Replicates

Number

1

☐ Vertical☒ Horizontal

Fill Layout

Delete

View: Identifiers

Direction

Multi Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	C1	STD1	S1	S2	S3	S4	S5	S6	S7	S8	S9	C2
B	P1	STD2	S10	S11	S12	S13	S14	S15	S16	S17	S18	N1
C	P2	STD3	S19	S20	S21	S22	S23	S24	S25	S26	S27	N2
D	P3	STD4	S28	S29	S30	S31	S32	S33	S34	S35	S36	N3
E	P4	STD5	S37	S38	S39	S40	S41	S42	S43	S44	S45	N4
F	P5	STD6	S46	S47	S48	S49	S50	S51	S52	S53	S54	N5
G	P6	STD7	S55	S56	S57	S58	S59	S60	S61	S62	S63	N6
H	C3	STD8	S64	S65	S66	S67	S68	S69	S70	S71	S72	C4

Cancel

&lt;

Back

Next

&gt;

## Method Selection

Available detection and preparation methods are displayed. To add detection or preparation methods, click-and-drag the method to the protocol, or select the method and click the Add button.

[Report Options](#)

Single Kinetic Area Scan Wavelength Scan

Select Method

- Shake
- Eject
- Load
- Pause
- ABS-MONO (s/n 1009)
  - 482nm
  - NewMethod 2
  - DNA 2 @977
  - DNA 2 @900
  - DNA 2 @320
  - DNA 2 @280
  - DNA 2 @260
  - Abs\_595nm
  - x\_Abs\_Mono\_495nm
  - 260nm
  - x\_Abs\_Mono\_495
  - x\_Abs\_Mono\_320nm
  - x\_Abs\_Mono\_382nm
  - x\_Abs\_Mono\_405nm
  - x\_Abs\_Mono\_280nm
  - x\_Abs\_Mono\_260nm

☒ Cartridge View

☒ Show only runnable Methods




Estimated Time

00:01:31

Group1

- x\_Abs\_Mono\_495
- Abs\_595nm



Cancel



Back

Next



## Data Reduction Page

Press F1 for more information about data reduction functions and formulas.

[Report Options](#)

General Settings  
Technique Type  
Labware Selection  
Layout Settings  
Method Selection  
**Data Reduction Page** >  
Concentration  
Cutoff  
Validation  
Output Settings

Group1

- A = x\_Abs\_Mono\_495
- B = Abs\_595nm

First Pass

[Add new Pass](#)

REDUCTION\_A1

Formula

Formula

B-A

Name of Data

Name of Units

495 nm

Notes

Apply Formula for Wells with Category

- ☒ Blank
- ☒ Control
- ☒ Empty
- ☒ Negative Control
- ☒ Positive Control
- ☒ Sample
- ☒ Standard

[Add new Item](#)

Cancel

<

Back

Next >

## Concentration

Choose a curve fitting method and configure response formulas and concentrations as desired.

### Standard Curve

☐ Use stored Standard Curve

☐ Copy standard curve definition from an other protocol to this protocol

Select

### Curve

Linear Regression

Number of standards 8

Extrapolation 0 %

### Y-Axis

Base REDUCTION\_A1 - (495 nm)

Type linear



### X-Axis

Name Concentration

Type linear



### Graph Setup

### Validation

Report Options

### Standards

	Response Formula	Concentration
1	STD1	1
2	STD2	2
3	STD3	3
4	STD4	4
5	STD5	5
6	STD6	6
7	STD7	7
8	STD8	8

Cancel

<

Back

Next >

## Cutoff

Cutoff configures qualitative evaluations that classify measured samples according to defined cutoff values. Up to ten groups of samples may be classified using cutoff formulas. Groups are separated by the cutoff formulas.

Basis of Evaluation

Report Options

Number of Groups

### Groups and Cutoff Formulas

Name		Max Value	Cutoff Formulas / Values
Group 2	<input type="text" value="High"/>		<input type="text" value="0.5"/>
Group 1	<input type="text" value="Low"/>		

Cancel

< Back

Next >

- General Settings
- Technique Type
- Labware Selection
- Layout Settings
- Method Selection
- Data Reduction Page
- Cutoff**
- Output Settings

## Cutoff

Cutoff configures qualitative evaluations that classify measured samples according to defined cutoff values. Up to ten groups of samples may be classified using cutoff formulas. Groups are separated by the cutoff formulas.

Basis of Evaluation

Report Options

Number of Groups

### Groups and Cutoff Formulas

	Name	Max Value	Cutoff Formulas / Values
Group 3	High		
		1.8	
Group 2	Normal		
		1.2	
Group 1	Low	0	

Cancel

<

Back

Next

>

## Validation

Configure up to 10 validation rules.

Number of Rules

4

Report Options

### Validation Rules

Basis of Evaluation	Rule
REDUCTION_A1 - (495 nm)	P1<0.8
If failed, show Text:	Failed, caused by Positive control 1
REDUCTION_A1 - (495 nm)	P4<0.8
If failed, show Text:	Failed, caused by Positive control 4
REDUCTION_A1 - (495 nm)	N1>0.5
If failed, show Text:	Failed, caused by Negative control 1
REDUCTION_A1 - (495 nm)	N4>0.5
If failed, show Text:	Failed, caused by Negative control 4

Cancel

<

Back





Next

>

## Output Settings

Select data output and printer options.

Perform after completing measurement(s)

- ☒  Export to Microsoft® Excel
- ☒  Show Result Viewer
- ☐  Create .XML and .dat data files
- ☐  Automatic Print out after measurement.

▷ Execute a program after protocol executes

▽ Save and run this protocol now

Save and run this protocol now.



Cancel



Back

Save

## Prepare to Run

Click Run to run the protocol.

Result Name

# Samples to measure

# Plates to read

Plate is lidded ☐

Eject Plate Carrier



Close Plate Carrier



Select Plate Orientation (click on Plate to view Layout)


☒ Landscape (Optimized)

☐ Portrait

☐ Opposite Landscape

☐ Opposite Portrait

Instrument



Run the selected protocol



### Optimization Features

☐ Do Labware Optimization before Run. Labware Optimization should be performed every Month to get best performance of your instrument. Last optimization was performed on 7/23/2007

Cancel

Run



## Run Protocol

Click Finish to view results and export data



Estimated Time 00:00:02

Elapsed Time 00:00:02

Current Method: Abs\_595nm

Raw Data Graph

Abs\_595nm  
Data

	1	2	3	4	5	6	7	8	9	10	11	12
A	2.071 OK	2.527 OK	3.419 OK	1.331 OK	0.301 OK	3.678 OK	2.923 OK	3.856 OK	2.343 OK	0.269 OK	3.477 OK	2.466 OK
B	2.278 OK	0.931 OK	0.327 OK	0.596 OK	1.308 OK	0.221 OK	0.116 OK	1.393 OK	2.176 OK	0.149 OK	0.755 OK	2.296 OK
C	3.896 OK	3.55 OK	0.443 OK	2.97 OK	0.55 OK	2.203 OK	3.945 OK	3.51 OK	2.099 OK	1.257 OK	0.728 OK	1.305 OK
D	3.195 OK	0.961 OK	2.277 OK	1.355 OK	1.94 OK	1.31 OK	2.227 OK	2.978 OK	3.044 OK	0.972 OK	3.04 OK	2.643 OK
E	1.667 OK	3.758 OK	3.681 OK	2.311 OK	0.27 OK	0.545 OK	3.191 OK	1.921 OK	1.771 OK	1.122 OK	1.435 OK	0.75 OK
F	3.234 OK	3.952 OK	3.305 OK	0.14 OK	0.324 OK	3.967 OK	0.367 OK	1.02 OK	0.203 OK	3.021 OK	1.4 OK	0.346 OK
G	1.943 OK	2.76 OK	3.452 OK	0.866 OK	1.922 OK	1.777 OK	3.251 OK	2.923 OK	0.51 OK	1.8 OK	1.302 OK	0.792 OK
H	2.521 OK	1.634 OK	3.239 OK	1.553 OK	2.066 OK	2.806 OK	3.534 OK	2.072 OK	3.526 OK	1.526 OK	2.121 OK	1.987 OK

Finish




## Data

View measurement results and graphs in the tabs below.

Raw data Edit Reduced Data GraphView

Abs\_595nm



Cancel

Next



Data

Concentration

Cutoff

Validation

## Data

View measurement results and graphs in the tabs below.

Raw data Edit Reduced Data GraphView

REDUCTION\_A1 - (595 nm)

	1	2	3	4	5	6	7	8	9	10	11	12
A	C1 2.071 OK	STD1 2.527 OK	S1 3.419 OK	S2 1.331 OK	S3 0.301 OK	S4 3.678 OK	S5 2.923 OK	S6 3.836 OK	S7 2.343 OK	S8 0.269 OK	S9 3.477 OK	C2 2.466 OK
B	P1 2.278 OK	STD2 0.931 OK	S10 0.327 OK	S11 0.596 OK	S12 1.308 OK	S13 0.221 OK	S14 0.116 OK	S15 1.393 OK	S16 2.176 OK	S17 0.149 OK	S18 0.755 OK	N1 2.296 OK
C	P2 3.896 OK	STD3 3.550 OK	S19 0.443 OK	S20 2.970 OK	S21 0.550 OK	S22 2.203 OK	S23 3.945 OK	S24 3.510 OK	S25 2.099 OK	S26 1.257 OK	S27 0.728 OK	N2 1.305 OK
D	P3 3.195 OK	STD4 0.961 OK	S28 2.277 OK	S29 1.355 OK	S30 1.940 OK	S31 1.310 OK	S32 2.227 OK	S33 2.978 OK	S34 3.044 OK	S35 0.972 OK	S36 3.040 OK	N3 2.643 OK
E	P4 1.667 OK	STD5 3.758 OK	S37 3.681 OK	S38 2.311 OK	S39 0.270 OK	S40 0.545 OK	S41 3.191 OK	S42 1.921 OK	S43 1.771 OK	S44 1.122 OK	S45 1.435 OK	N4 0.750 OK
F	P5 3.234 OK	STD6 3.952 OK	S46 3.305 OK	S47 0.140 OK	S48 0.324 OK	S49 3.967 OK	S50 0.367 OK	S51 1.020 OK	S52 0.203 OK	S53 3.021 OK	S54 1.400 OK	N5 0.346 OK
G	P6 1.943 OK	STD7 2.760 OK	S55 3.452 OK	S56 0.866 OK	S57 1.922 OK	S58 1.777 OK	S59 3.251 OK	S60 2.923 OK	S61 0.510 OK	S62 1.800 OK	S63 1.302 OK	N6 0.792 OK
H	C3 2.521 OK	STD8 1.634 OK	S64 3.239 OK	S65 1.553 OK	S66 2.066 OK	S67 2.806 OK	S68 3.534 OK	S69 2.072 OK	S70 3.526 OK	S71 1.526 OK	S72 2.121 OK	C4 1.987 OK

Cancel

Next



Data  
Concentration  
Cutoff  
Validation

## Cutoff

View cutoff groups or edit the cutoff formulas if reevaluating the results is desired.

Layout View List View Parameters Edit

### Results

☐ Show Well Identifier

☒ Show Status

	1	2	3	4	5	6	7	8	9	10	11	12
A	High OK	High OK	High OK	High OK	Low OK	High OK	High OK	High OK	High OK	Low OK	High OK	High OK
B	High OK	High OK	Low OK	High OK	High OK	Low OK	Low OK	High OK	High OK	Low OK	High OK	High OK
C	High OK	High OK	Low OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK
D	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK
E	High OK	High OK	High OK	High OK	Low OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK
F	High OK	High OK	High OK	Low OK	Low OK	High OK	Low OK	High OK	Low OK	High OK	High OK	Low OK
G	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK
H	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK	High OK

Cancel

< Back

Next >



Data  
Concentration  
Cutoff  
Validation >

## Validation

View results for each validation rule or edit the rules if reevaluating the results is desired.

Parameters Edit

Result List

Rule Number	Result	FailText	Basis of Evaluation	Rule
1	Pass		REDUCTION_A1 - (595 nm)	C1>2
2	Pass		REDUCTION_A1 - (595 nm)	P1>1.5
3	Pass		REDUCTION_A1 - (595 nm)	N1>1.5

Cancel

< Back

Close



Data  
Concentration  
Cutoff  
Validation >

## Validation

View results for each validation rule or edit the rules if reevaluating the results is desired.

Parameters

Edit

Number of Rules 3

Report Options

### Validation Rules

Basis of Evaluation	Rule
REDUCTION_A1 - (595 nm)	C1>2
If failed, show Text:	Caused by Control 1 failure
REDUCTION_A1 - (595 nm)	P1>3
If failed, show Text:	Caused by P1 failure
REDUCTION_A1 - (595 nm)	N1>1.5
If failed, show Text:	Caused by N1 failure

Cancel

< Back

Close



Data  
Concentration  
Cutoff  
Validation >

## Validation

View results for each validation rule or edit the rules if reevaluating the results is desired.

Parameters

Edit

Result List

Error Validation failed.

Rule Number	Result	FailText	Basis of Evaluation	Rule
1	Pass		REDUCTION_A1 - (595 nm)	C1>2
2	Fail	Caused by P1 failure	REDUCTION_A1 - (595 nm)	P1>3
3	Pass		REDUCTION_A1 - (595 nm)	N1>1.5

Cancel

< Back

Close

Multimode Analysis Software

FileActionsHelp

Delete AllDeleteViewPrint

Protocols

Detection Methods

Results

Labware

Instruments

Trash

Results Selection List

Result Name	Measured	Evaluated	Protocol Name	Application Type
20070816-173	2007/8/16 下午 05:34:13	2007/8/16 下午 05:43:23	Validation 595nm	Analysis
validation 595n	2007/8/15 下午 04:16:59	2007/8/15 下午 04:17:02	Validation 595nm	Analysis
20070815-160	2007/8/15 下午 04:07:41	2007/8/15 下午 04:09:57	NewProtocol 7	Analysis
20070815-152	2007/8/15 下午 03:28:35	2007/8/15 下午 03:28:37	Abs_595nm	Analysis
20070815-141	2007/8/15 下午 02:17:28	2007/8/15 下午 02:27:34	NewProtocol 6	Analysis
20070815-132	2007/8/15 下午 01:28:29	2007/8/15 下午 01:28:31	Abs_595nm	Analysis
20070815-125	2007/8/15 下午 12:57:20	2007/8/15 下午 12:57:22	Abs_595nm	Analysis
20070814-153	2007/8/14 下午 03:31:57	2007/8/14 下午 03:32:02	NewProtocol 5	Analysis
20070814-152	2007/8/14 下午 03:23:07	2007/8/14 下午 03:23:40	Abs_580	Analysis
20070810-164	2007/8/10 下午 04:49:56	2007/8/10 下午 04:51:07	x_Abs_Gen_96well	Analysis
20070810-164	2007/8/10 下午 04:48:36	2007/8/10 下午 04:48:51	x_Abs_Gen_96well	Analysis
20070725-224	2007/7/25 下午 10:44:07	2007/7/25 下午 10:44:11	Lum_1000	Analysis
20070725-223	2007/7/25 下午 10:36:37	2007/7/25 下午 10:36:41	x_Multi_Luminescence 96well 2.5min	Analysis
20070607-163	2007/6/7 下午 04:37:17	2007/6/7 下午 04:37:24	NewProtocol 4	Analysis
20070607-161	2007/6/7 下午 04:16:17	2007/6/7 下午 04:16:23	x_Multi_Rhodamine Int Bottom 96well 2min 5x5pts	Analysis
20070605-180	2007/6/5 下午 06:05:10	2007/6/5 下午 06:05:43	NewProtocol 3	Analysis
20070605-175	2007/6/5 下午 05:56:56	2007/6/5 下午 05:57:00	NewProtocol 3	Analysis
20070605-133	2007/6/5 下午 01:40:36	2007/6/5 下午 01:40:42	NewProtocol 2	Analysis
20070605-131	2007/6/5 下午 01:19:59	2007/6/5 下午 01:20:05	NewProtocol 2	Analysis
20070605-115	2007/6/5 上午 11:57:48	2007/6/5 上午 11:57:55	x_Multi_Europium Chelate 384well 4min	Analysis

Result NameLike

DateBetween2007/ 8/14and2007/ 8/16

Protocol nameLike

GoClear

LTTW01A0006

PARADIGM - Simulated

TonBottom

Jack YE...ParadigmParadig...Multim...Validate...Validate...CE100%

下午 05:52



謝謝