

WELGENE 威健生技

全方位基因體及轉錄體分析

李彥樑
行銷暨業務經理

台北市南港軟體園區 F 棟 12 樓
+886-2-6616001 ext. 223
jacklee64@welgene.com.tw



Welgene Biotech. Co. Ltd
<http://www.welgene.com.tw>



威健生技 高通量服務的領導者



Office Area



Service Lab

Biotech Building, Nan-Kang Business Park



Quality is Assured

Welgene Biotech. Co. Ltd
<http://www.welgene.com.tw>

Agilent Certification

- Agilent certified service provider since 2007

Agilent Technologies | Genomics | Login | Contact Us

Certified Service Providers - Asia Pacific

CSP	Certified Applications Provided							
	1-color Gene Expression	2-color Gene Expression	microRNA	ChIP-on-chip	CGH with enzymatic labeling	CGH with non-enzymatic or ULS	CGH with Automated High-throughput	Target Enrichment
CERI	+	+	+	+	+	+		
The Chinese University of Hong Kong	+	+	+	+	+	+		
Welgene Taipei, Taiwan	+	+	+	+	+	+		
Genotypic Technology Pte Ltd.	+	+	+	+	+			
Hokkaido System Science Co	+	+	+	+	+	+		
National Yang Ming University	+	+						
Welgene	+	+	+	+	+	+		

Welgene Biotech has implemented the complete Agilent microarray platform to ensure the best array experiment quality. The array service quality has certified by ISO17025. The key steps in the workflow include: Consult for array experiment design, Capture cell by LCM, Extract DNA/RNA for customer, Performs meDIP for customer, Design Custom GE/CGH/CIP-on-chip array design, Analyze data with GeneSpring GX and DNA Analytics. Services offered to customers in Taiwan and Hong Kong

ISO17025 Certification

- ISO certification since 2006 (SBIR support)

認證依據：ISO/IEC 17025：2005

初次認證日期：九十五年八月二十九日

認證範圍：測試領域，如續頁

實驗室主管：林怡杏

認可事項

14.01 生物科技

核糖核酸,組織,細胞株

B999 基因表現

RNA Isolation (Sr-03-A)

Low RNA Input Fluorescent Linear Amplification & Agilent Oligo Microarray Hybridization Spike-In Control(Sa-03-A)

測試組與對照組做相對量測試

報告簽署人：林怡杏,林燦濤,黃信智

- **Hepatogastroenterology. 2006 Jul-Aug;53(70):484-90.** Distinct gene expression profiles in gastric epithelial cells induced by different clinical isolates of Helicobacter pylori--implication of bacteria and host interaction in gastric carcinogenesis.
- **Ann Surg Oncol. 2006 Nov;13(11):1474-84.** A gene expression profile for vascular invasion can predict the recurrence after resection of hepatocellular carcinoma: a microarray approach.
- **Fertility and Sterility 2006, Vol. 86(6):1650-1658,** Identification of ten novel genes involved in human spermatogenesis by microarray analysis of testicular tissue.,
- **Oncol Rep. 2006 Apr;15(4):919-26.** A homologue of the Drosophila headcase protein is a novel tumor marker for early-stage colorectal cancer.
- **Oncogene. 2007.** Molecular Signatures of Metaplastic Carcinoma of Breast by Large-Scale Transcriptional Profiling: Identification of Genes Potentially Related to Epithelial-Mesenchymal Transition..
- **BMC Genomics. 2008 Feb 29;9:109** Altered expression patterns of lipid metabolism genes in an animal model of HCV core-related, nonobese, modest hepatic steatosis.
- **Breast Cancer: Basic and Clinical Research. 2008.** Comparison and Identification of Estrogen-Receptor Related Gene Expression Profiles in Breast Cancer of Different Ethnic Origins.

2009 Labeled microRNA pull-down assay system: an experimental approach for high-throughput identification of microRNA-target mRNAs

Nucleic Acids Research, 2009, Vol. 37, No. 10 e77

2009 Statistical identification of gene association by CID in application of constructing ER regulatory network

BMC Bioinformatics 2009, 10:85 doi:10.1186/1471-2105-10-85

2009 Interactions between octaarginine and U-937 human macrophages: global gene expression profiling, superoxide anion content, and cytokine production.

J Control Release. 2009 Nov 3;139(3):197-204. Epub 2009 Jul 18.




2009 Hepatic inflammation mediated by hepatitis C virus core protein is ameliorated by blocking complement activation.

BMC Med Genomics. 2009 Aug 8;2:51.

2010 A time-course study of gene responses of chicken granulosa cells to Salmonella Enteritidis infection.


Vet Microbiol. 2010 Jan 25. [Epub ahead of print]

服務內容



分析技術

Microarray




Agilent Technologies


Q-RT-PCR Array



OpenArray™ architecture



NGS



Illumina Solexa




AB SOLiD



Roche 454



Agilent Technologies



全方位的核酸分析服務

SNP Analysis

(Q-PCR array / Transcriptome-Seq)

aCGH / CNV

(microarray)

DNA Methylation

(microarray / MeDIP-Seq)

Gene Expression

(microarray / DGE / RNA-Seq)

Alternative Splicing

(Transcriptome-Seq / Deep-Seq)


miRNA Expression

(microarray / miRNA-Seq)


Transcription Factor Binding

Histone Modification

(ChIP-on-chip / ChIP-Seq)



Integration of Multi-omics Research




Applications

- ✚ Genomics Study (SNP, CNV)
- ✚ Transcriptomics Study (mRNA / miRNA expression, Alternative Splicing)
- ✚ Transcription Factor Binding (ChIP assay)
- ✚ Epigenomics Study (Methylation, Histone Modification)


Genomics Study

Normal diploid genome



Copy Number Variation (CNV):

- Copy number variation (CNV): DNA segments in which copy-number varies between two or more genomes
- Ranges from a thousand to millions of DNA bases in size
- CNVs have been associated with susceptibility to disease, complex behavioral traits, and other phenotypic variability
- Identifying statistically significant CNVs is important in understanding the underlying mechanism of disease and disease susceptibility




Single Nucleotide Polymorphisms (SNPs):

- DNA sequence variations that occur when a single nucleotide in genome is altered
- Variation must occur in at least 1% of the population to be considered a SNP
- Occur every 100 to 300 bases

aCGH / CNV Array分析

Array CGH




探針


染色體

尋找癌症或疾病相關的染色體變異點

Human colon carcinoma cell line HT29



PWS/AS deletion



Database of Genomic Variants

Hosted by:
The Centre for



A curated catalogue of structural variation in the human genome

[About The Project](#) | [Genome Browser](#) | [Download](#) | [Lit](#)

Please select genome assembly:

View Data by Chromosome

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 X

Keyword Search

Search

Exact Match? Yes No

Examples: clone name, accession number, cytoband or gene

Genet Med 2008;10(6):415–429.

article

Enhanced detection of clinically relevant genomic imbalances using a targeted plus whole genome oligonucleotide microarray

Erin L. Baldwin, PhD, Ji-Yun Lee, PhD, Douglas M. Blake, BS, Brian P. Bunke, BS, Chad R. Alexander, BS, Amy L. Kogan, BS, David H. Ledbetter, PhD, and Christa L. Martin, PhD


wellcome trust
sanger
institute

RSS

Human Genetics | Model Organisms | Pathogens | Bioinformatics | Sequencing

The Copy Number Variation (CNV) Project

Genetic diseases are caused by a variety of of large chunks of DNA sequence consisting mutation has often been overlooked in previc disease is caused by copy number variation occur in families result from these kinds of r malaria. The contribution of CNV to the com Mutations (of any type) that increase the ch developing a common disease are also likel frequencies in apparently healthy population...



High-resolution array CGH analysis of salivary gland tumors reveals fusion and amplification of the FGFR1 and PLAG1 genes in ring chromosomes.


Persson F, Winnes M, Andrén Y, Wedell B, Dahlenfors R, Asp J, Mark J, Enlund F, Stenman G.

Lundberg Laboratory for Cancer Research, Department of Pathology, Sahlgrenska University Hospital, Göteborg University, Göteborg, Sweden.

Welgene Biotech. Co. Ltd
<http://www.welgene.com.tw>

Requirement for Specimens

- DNA
 - Quantity
 - ≥ 4 ug
 - Quality
 - $OD_{260/280} > 1.7$,
 - size > 4 kb




aCGH /CNV Result

Analysis Settings

Aberration Algorithm : ADM-2
 Threshold : 6.0
 Centralization : ON
 Bin Size : 10
 Centralization Threshold : 6.0
 Fuzzy Zero : ON
 Combine Replicates (Intra Array) : OFF
 Genome : hg18
 Aberration Filters : minProbes = 3 AND minAvgAbsLogRatio = 0.4 AND maxAberrations = 100 AND percentPenetrance = 0
 Feature Level Filters : NONE
 Expand Non Unique Probes : OFF
 Genomic Boundaries : Not Applied

Genome Overview




Aberration	CytoBand	ChrName	ProbeName	Start	Stop	Descriptor	Genes	hg18_CNVHs	hg18_(Human hg18 cytoband)	LogRatio
1.1	p32.1	chr1	A_16_P00	59012360	59012417	Unknown			1p32.1	0.486557
1.2	p32.1	chr1	A_14_P20	59018512	59018571	Unknown			1p32.1	0.647439
1.3	p32.1	chr1	A_14_P12	59019101	59019160	Homo sapi	JUN		1p32.1	0.487863
1.4	p32.1	chr1	A_14_P13	59021301	59021350	Homo sapi	JUN		1p32.1	0.911086
1.5	p32.1	chr1	A_14_P10	59022152	59022210	Homo sapi	JUN		1p32.1	0.757187
1.6	q21.1	chr1	A_16_P00	1.47E+08	1.47E+08	Homo sapiens neurobl	CNV_3315	CNV_8288	1q21.1	-0.50044
1.7	q21.1	chr1	A_16_P00	1.47E+08	1.47E+08	Unknown	CNV_3315	CNV_8288	1q21.1	-0.77127
1.8	q21.1	chr1	A_16_P00	1.47E+08	1.47E+08	Unknown	CNV_3315	CNV_8288	1q21.1	-0.69443
1.9	q21.1	chr1	A_14_P13	1.47E+08	1.47E+08	Unknown	CNV_3315	CNV_8288	1q21.1	-0.32059
1.1	q21.1	chr1	A_16_P35	1.47E+08	1.47E+08	Unknown	CNV_3315	CNV_8288	1q21.1	-1.26298
1.11	q22	chr1	A_16_P00	1.53E+08	1.53E+08	Homo sapiens glucosic	CNV_3315	CNV_3315	1q22	-0.51631
1.12	q22	chr1	A_16_P35	1.53E+08	1.53E+08	Homo sapiens glucosic	CNV_3315	CNV_3315	1q22	-0.12569
1.13	q22	chr1	A_16_P00	1.53E+08	1.53E+08	Homo sapiens glucosic	CNV_3315	CNV_3315	1q22	-0.45874
1.14	q22	chr1	A_16_P00	1.53E+08	1.53E+08	Homo sapiens glucosic	CNV_3315	CNV_3315	1q22	-0.38016
1.15	q22	chr1	A_16_P15	1.53E+08	1.53E+08	Homo sapiens glucosic	CNV_3315	CNV_3315	1q22	-0.80982

Technician: _____ Date: __/__/__ Supervisor: _____ Date: __/__/__


Text Summary Report for Sample GC-05_US85003611_251469359030_S01

Event No	Chr	Cytoband	#Probes	Amp/Del	P-value	Annotations
1	chr1:59012160-59022352	p32.1	5	0.658026	3.63E-10	JUN, 1p32.1
2	chr1:147386718-147478054	q21.1	5	-0.677868	3.42E-10	CNV_3313, CNV_2326, CNV_1328...
3	chr1:153451703-153463890	q22	5	-0.458142	1.07E-10	CNV_3315, CNV_4257, CNV_3315...




SNP Q-PCR Array System

基因分型儀器平臺 (TaqMan[®] assay)



OpenArray[™] architecture




Service Flow

Welgene Service Lab



1. Search Taqman assay using SNP ID

2. DNA samples



Order SNP Array assay
Turn around time: 6 weeks

DNA QC


Samples mixed with Master Mix (containing ROX), and loaded to OpenArray

PCR

Data readout


Data QC

Assay Format	Total sample number
16 SNP	1,269
32 SNP	837
64 SNP	874
128 SNP	836




Requirements for input samples

- DNA
 - In nuclease-free water, NOT in T.E
 - Quantity:
 - 400 ng (w/t WGA)
 - Concentration: ~50ng/ul
 - Quality:
 - $OD\ 260/280 > 1.7$
 - $OD\ 230/280 > 1.7$
 - Agarose gel electrophoresis
 - $> 10\ kb$




Assay Result

	A	B	C	D	E	F	G
1	Open Array, Serial Number	Sample, SampleID	Sample, Description	C__29086771_20, Genotype	C__3168989_10, Genotype	C__8376164_10, Genotype	C__3197300_10, Genotype
2	CAD83	NA17102		No Call	G G	A A	G A
3	CAD83	NA17103		No Call	G A	G A	G G
4	CAD83	NA17104		G G	G G	No Call	A A
5	CAD83	NA17105		G G	G A	A A	G A
6	CAD83	NA17106		No Call	No Call	A A	G A
7	CAD83	NA17107		G A	G A	G A	A A
8	CAD83	NA17108		G G	G A	G A	G A
9	CAD83	NA17109		G A	G G	G A	G A
10	CAD83	NA17111		G A	G A	G A	G A
11	CAD83	NA17115		G G	G A	A A	G A
12	CAD83	NA17117		G A	G G	G A	G A
13	CAD83	NA17125		G G	G A	G G	G A




62 assays, average call = 97.8%




Data points correspond to the two homozygote genotypes (blue, red), the heterozygote genotype (green), the no template control (NTC) (black) for the unamplified DNA samples, and the WGA amplified samples (orange). Samples that fall within the gray box demarcated by the ⊗ symbols are not called. In all cases the call rate was 99% and accuracy was equal to 100%. a.u., arbitrary units.

Transcriptomics Study




Gene Expression

Gene expression is the process by which information from a gene is used in the synthesis of a functional gene product.



miRNA expression:

- ✚ 20~24 nt short single strand RNA molecule
- ✚ Endogenous non-coding sequence
- ✚ Acting as translation inhibitor by degrading mRNA or blocking protein synthesis



Alternative Splicing:

Alternative splicing occurs as a normal phenomenon in eukaryotes, where it greatly increases the diversity of proteins that can be encoded by the genome; in humans, over 80% of genes are alternatively spliced. There are numerous modes of alternative splicing observed, of which the most common is exon skipping. In this mode, a particular exon may be included in mRNAs under some conditions or in particular tissues, and omitted from the mRNA in others.



Agilent Technologies

60-mer Oligo Expression Array



Human Whole Genome Array 40,000 genes and transcripts



Mouse Whole Genome Array 41,000 genes and transcripts



Rat Whole Genome Array 41,000 genes and transcripts

另有其它物种

Monkey

Pig

Dog

Cat

Arabidopsis

Rice

Yeast


Zebrafish

M. grisea

C. elegans

Xenopus laevis

More...



Welgene Biotech. Co. Ltd
<http://www.welgene.com.tw>

市面最新版的 miRNA array (miRBase V14)




Figure 1. Components of the Agilent miRNA microarray probe design. An unmodified microarray probe (black) is a synthesized sequence that hybridizes to the target miRNA (red). Probes are anchored to the glass slide surface by a stilt (brown). **A.** Inclusion of a G residue (black) to the 5' end of the hybridization sequence complements the 3' end C residue (yellow) introduced in labeling. This additional G-C pair in the probe-target interaction region stabilizes targeted miRNAs relative to homologous RNAs. Additionally, all probes contain a 5' hairpin (blue), abutting the probe-target region, to increase target and size miRNA specificity. **B.** Destabilization of probes that are too stable. For probes requiring it, reduction of probe-target base-pairing is achieved through sequential elimination of base pairing from the 5' end of the miRNA.

Human rel12 miRNA Array

955 human-relative miRNA

Mouse rel14 miRNA Array



690 mouse-relative miRNA

Rat rel14 miRNA Array

388 rat-relative miRNA


Dual-Mode Array System

Two color



Red/Green

One color



Green/Green

Fold change of gene expression

Service Procedures

1. RNA QC Test

-check by Agilent 2100 Bioanalyzer

2. Sample Amplification & Labeling







-process with full sets of Agilent kit

3. Array Hybridization, Wash, & Scanning


-Agilent scanner and software provide superior resolution and extended signal range

4. Data Output

-easy to use data presented by Welgene



Array Report



GeneName	SystematicName	Description	PValueLogRatio	Normalized Exprssion Ratio
SPINK1	NM_003122	Homo sapiens serine peptidase inhibitor, Kazal type 1 (SPINK1), mRNA [NM_003122]	0.000	536.450
RAB3B	ENST00000371655	Ras-related protein Rab-3B. [Source:Uniprot/SWISSPROT;Acc:P20337] [ENST00000371655]	0.000	155.765
KCNC2	NM_139136	Homo sapiens potassium voltage-gated channel, Shaw-related subfamily, member 2 (KCNC2), transcript variant 1, mRNA [NM_139136]	0.000	107.736
GUCY2C	NM_004963	Homo sapiens guanylate cyclase 2C (heat stable enterotoxin receptor) (GUCY2C), mRNA [NM_004963]	0.000	75.554
MUC13	NM_033049	Homo sapiens mucin 13, cell surface associated (MUC13), mRNA [NM_033049]	0.000	73.358
LOC344887	BX640843	Homo sapiens mRNA; cDNA DKFZp686B14224 (from clone DKFZp686B14224). [BX640843]	0.000	67.974
REG3A	NM_138938	Homo sapiens regenerating islet-derived 3 alpha (REG3A), transcript variant 2, mRNA [NM_138938]	0.000	67.390

qRT-PCR Service (for RNA or DNA)





LightCycler



Mx-3000p

TaqMan Assay and SyBR Assay are both available

RNA-Seq



Illumina Solexa





AB SOLiD




Roche 454

Alternative Splicing / Point Mutation Analysis

每個 RNA 樣品需要準備 20 ug, RNA品質需 RIN >7
 進行1-2 G bases的定序量
 分析方式需特別討論




Epigenomics Study



DNA Methylation:

- ✚ The CG island is a short stretch of DNA in which the frequency of the CG sequence is higher than other regions. It is also called the CpG island, where "p" simply indicates that "C" and "G" are connected by a phosphodiester bond.
- ✚ CpG islands are often located around the promoters of housekeeping genes (which are essential for general cell functions) or other genes frequently expressed in a cell. At these locations, the CG sequence is not methylated. By contrast, the CG sequences in inactive genes are usually methylated to suppress their expression.

Histone Modification:




Histones are subject to a wide variety of posttranslational modifications including but not limited to, lysine acetylation, lysine and arginine methylation, serine and threonine phosphorylation, and lysine ubiquitination and sumoylation (Vasquero 2003). The existence of these modifications and recognition modules led to a well established “histone code” hypothesis proposed by Strahl and Allis (2000). Overall, posttranslational modifications of histones create an epigenetic mechanism for the regulation of a variety of normal and disease-related processes.


Chromatin ImmunoPrecipitation

Genome wide location analysis to examine three different types of targets:

- i) Transcription factors
- ii) Cofactors
- iii) Histone modifiers
- iv) DNA methylation



Antibody:
specific to TF or methylated DNA




***In vivo* protein-DNA complexes**

IP chromatin to capture crosslinked protein-DNA complexes
 Sonicate complexes to fragment DNA
 Reverse and remove crosslinks

Purify and amplify DNA fragments
 =
 ENRICHED POOL of protein-bound DNA

ChIP-on-chip /ChIP-Seq Concept

Regulatory proteins bind to promoter DNA regions *in vivo*




Hybridize to microarray for detection




protein-bound DNA fragments

ChIP Analytics output in UCSC browser view



1. Cross-link protein-DNA complexes
2. Lyse cells and sonicate DNA
3. IP chromatin to capture and purify bound DNA
4. Release and amplify DNA fragments

Sequencing IP-DNA with NGS



Promoter Array Design


Goals:

- Construct a database of high-quality probes spanning the genome.
- Provide tools to select probes onto arrays for particular applications.

Probe optimization criteria:

- **Uniqueness (low homology)**
- **T_m**
- **Self-structure**
- **Balances spacing with performance**

7.5 kb promoter region— 5.5 kb upstream, 2 kb downstream



60-mer oligonucleotide probes—critical for the sensitivity and specificity

Average probe spacing parameters have been specifically optimized

Repeat regions are masked to significantly reduce non-specific noise.


Methylation (CpG Island) Array Design




- Tiled arrays with 60-oligomer probes spaced ~100 bp apart for 28,500 islands in human and 16,030 in mouse

- Uses the CpG island probes defined by Gardiner-Garden and Frommer from UCSC Hg18

- Compatible with methylated DNA immunoprecipitation method (Keshet et. al, 2006) and possibly other methods




客戶端製備 chromatin-immunoprecipitation-DNA



需準備 5×10^7 至 1×10^8 細胞 · 以 formaldehyde 處理細胞以 cross-link 蛋白質及核酸



磁珠與 ChIP 抗體結合



溶解細胞，並以超音波將 chromatin 片段化至長度 100bp-500bp

進行 chromatin-immunoprecipitation 實驗 · 解除 cross-link，並純化 DNA

以基因特異引子確認 ChIP 實驗成功


DNA 樣品需求

ChIP DNA, ~500 ng, 長度分布在 100~500 bp

Input DNA, ~ 8 μ g, 長度分布在 100~500 bp

ChIP 基因特異引子, 5 μ m, 20 μ l

ChIP-DNA、Input-DNA 及基因特異引子送交威健實驗室




確認 DNA 長度分布




以基因特異引子進行定量 PCR 確認 ChIP 實驗成功



補平 DNA 兩端，並連接上 linker·進行 LM-PCR 以增幅 ChIP- 及 Input-DNA 量



反應，螢光標幟 ChIP- 及 Input-DNA



進行 hybridization 及 wash 步驟




輸出，並分析 binding element


Slide Name: US22502707_251390210429_S01

Number of Probes	44290
Significant Probes	426
Bound Probes	440
Derivative Log ratio spread	0.2474
Normalized Red/Green Correlation (R ²)	0.4031

IP vs WCE





Enrichment ratio (logscale) vs Intensity



U	V	W	X	Y
P[Xbar]	Is In BoundRegion	Exclude	Primary An	Secondary
1	0	0	HsCGHBri	NONE
1	0	0	NegativeCo	NONE
0.739408	0	0	FAM79A	ref:FAM79
0.793869	0	0	DENND1B	ref:DENND1
0.001788	0	0	TOR1AIP2	ref:TOR1A
0.880848	0	0	PDE6B	ref:PDE6B:
0.363265	0	0	DKFZP434	ref:DKFZP:
1	0	0	LACC:GD2	NONE
0.759735	0	0	MTCH1	ref:MTCH1
0.565234	0	0	ESR1	ref:ESR1:-5
0.085293	0	0	MRPS24	ref:MRPS2:
0.706946	0	0	FLJ35779	ref:FLJ357
0.133019	0	0	FRS3-C6or	ref:FRS3:-1
0.356215	0	0	DHFRL1	ref:DHFRL
0.089522	0	0	RQCD1	ref:RQCD1
0.017177	0	0	PTGS2	ref:PTGS2:

ChIP-on-chip Result



Polycomb (Suz12) Covers Promoters and Genes

■	Pol II
■	Suz12

Methylation Array Result






Figure 3. Methylation comparison of the HoxA gene cluster between tumor and normal lung cell lines. Methylated DNA was immunoprecipitated and labeled from genomic DNA (as previously described) from the small cell lung carcinoma NCI-H69 cell line (**Figure 3A**) and the normal embryonic lung fibroblast HEL299 cell line (**Figure 3B**). For each line, samples of both Cy5-labeled mDIP-enriched DNA and Cy3-labeled input genomic DNA were competitively hybridized to an Agilent Human CpG Island Microarray (P/N G4492A). Microarrays were processed, scanned on an Agilent Scanner (P/N G2565BA), and data normalized to a median log₂ ratio of zero. The methylation profile for a portion of the Hox locus on Chromosome 7 is shown. Every point represents an individual CpG island microarray probe. Each colored line (orange for tumor, blue for normal) represents a three-point moving average of the log₂ ratio. Regions of more methylation (red) and less methylation (green) are indicated. Highly methylated regions have log₂ ratios significantly above zero while less methylated regions have log ratios significantly below zero. Profiles clearly show significant differences between the normal and cancer samples at HoxA6, HoxA7, and HoxA9 promoter regions.

Histone Modification Result




- Histone acetylation enriched at promoter regions and transcriptional start sites of active genes.
 - Acetylation at Gcn5 targets H3K9 and H3K14 and Esa1 targets H4 at Lys 5,8,12 and16.
 - Model: transcriptional activators recruit Gcn5 and Esa1 to promoters of genes upon their Activation (Robert et al., 2004) .

ChIP-Seq




Based on data: from Barski et al. 2007, Cell 129, 823-837.

Integration of Multi-omics Research



進階數據分析服務




Find Significantly Expressed Genes


ProbeName	GeneSymbol	Common name	Description
A_23_P64873	DCN	NM_001920	Homo sapiens decorin (DCN), transcript variant A1, mRNA [NM_001920]
A_23_P421401	PDGFRB	NM_002609	Homo sapiens platelet-derived growth factor receptor, beta polypeptide (PDGFRB), mRNA [NM_002609]
A_23_P209625	CYP1B1	NM_000104	Homo sapiens cytochrome P450, family 1, subfamily B, polypeptide 1 (CYP1B1), mRNA [NM_000104]
A_23_P211233	COL6A2	NM_001849	Homo sapiens collagen, type VI, alpha 2 (COL6A2), transcript variant 2C2, mRNA [NM_001849]
A_23_P164057	MFAP4	NM_002404	Homo sapiens microfibrillar-associated protein 4 (MFAP4), mRNA [NM_002404]
A_23_P142533	COL3A1	NM_000090	Homo sapiens collagen, type III, alpha 1 (Ehlers-Danlos syndrome type IV, autosomal dominant) (COL3A1), mRNA [NM_000090]
A_23_P0063	LUM	NM_002345	Homo sapiens lumican (LUM), mRNA [NM_002345]
A_23_P0068	CPM	NM_001874	Homo sapiens carboxypeptidase M (CPM), transcript variant 1, mRNA [NM_001874]
A_23_P0068	SPON2	NM_012445	Homo sapiens spondin 2, extracellular matrix protein (SPON2), mRNA [NM_012445]
A_23_P0068	PRRX2	NM_016307	Homo sapiens paired related homeobox 2 (PRRX2), mRNA [NM_016307]
A_23_P0068	COL3A1	NM_000090	Homo sapiens collagen, type III, alpha 1 (Ehlers-Danlos syndrome type IV, autosomal dominant) (COL3A1), mRNA [NM_000090]
A_23_P0068	GPM6B	NM_001001996	Homo sapiens glycoprotein M6B (GPM6B), transcript variant 2, mRNA [NM_001001996]
A_23_P0033	PDGFRA	NM_006206	Homo sapiens platelet-derived growth factor receptor, alpha polypeptide (PDGFRA), mRNA [NM_006206]

Network Analysis



GO Analysis




BioPax, NLP, MeSH



整合性數據分析



Find overlapped genes



Build possible gene interaction network



**Thank you for your attending
Wish you have a day**



Welgene Biotech Co., Ltd.

Welgene Biotech. Co. Ltd
<http://www.welgene.com.tw>