




全自動毛細管DNA遺傳分析儀 - ABI PRISM 3100 的原理及應用

美商應用生命系統(股)公司
Field Application Specialist
徐英誠

Terminator Sequencing Sample Prep




1 reaction within which each terminator is labeled with a different dye



Sample is purified by ethanol precipitation
to remove unincorporated ddNTP's before loading

Sanger Dideoxy Cycle Sequencing


Performed in a thermocycler with a thermostable Taq polymerase



AB Applied Biosystems

Dye Terminator Cycle Sequencing Step

(1) Template +Primer +Sequencing reagents (with dye terminator)



Advantages of Capillary Electrophoresis



- More efficient heat dissipation than slab gels
 - ◆ higher run voltages
 - ◆ faster run times
 - ◆ **higher sample throughput**
- Automation
 - ◆ elimination of manual operations
 - ◆ increased run to run consistency and reliability
- Sensitivity
 - ◆ less DNA per sample
 - ◆ but more strict sample requirement(ex: aware of salt contamination of the seq. product....., etc)



ABI PRISM® 3100 Overview




- **16 Capillaries**
- Automated polymer filling, sample injection, detection & analysis
- **Supports Sequencing & Fragment Analysis**
- Supports capillary lengths 36 cm and 50 cm
 - ◆ ^{New} 22 cm for SNP genotyping and 80 cm for long sequence application (up to 1000 bp)
- **Active temperature control**
 - ◆ Capable of subambient cooling (18°C to 65 °C)
- Autosampler supports 96 & 384
- **Up to 24 hours unattended operation**
- PC running Windows XP




Process of Typical Run

AB Applied Biosystems


1. Sample Preparation




3. Beginning the Run




2. Software Setup




4. Electrophoresis




5. Excitation and Detection



7. Data Processing: during running



6. Data Collection




8. Automatic Data Extraction and Data Analysis




AB Applied Biosystems

9. Viewing the Data: on the GeneMapper software or Sequencing Analysis software




ABI PRISM® 3100 : Doors Open


AB Applied Biosystems



Autosampler Close-up



- Supports up to two 96- and or 384-well micrometer plates
- Automated sample injection
- Autosampler will stay up if power is lost




Instrument Open- Oven and Array



- Peltier heat pump used to heat or cool oven
 - * can cool down quickly
 - * facilitates subambient cooling
- Force air oven with two internal fans to provide superior thermal control and uniformity
- Array combs that clip into oven to provide reproducible array spatial alignment
- No more tape!!

3130 System Optics




- Laser source:
 - * 25 mW single argon ion laser
 - * running at 15 mW
 - * simultaneous dual side illumination

- Detection:
 - * In - capillary detection
 - * two- dimensional (spectral & spatial dimensions) data collected by the CCD camera



Schematic of “In-capillary” Detection




Close-up image of capillary array detection cell with a stylized schematic representation of in-capillary detection.



Dual Side Illumination

AB Applied Biosystems



AB Applied Biosystems

Electrokinetic Injection

AB Applied Biosystems

- Capillary and electrode are placed into the sample
- Voltage is applied
- “-” charged DNA enters the capillary as it migrates toward the “+” electrode at the other end of the capillary

Because of the nature of electrokinetic injection, sample preparation is critical to the data quality. Ex: salt contamination, template/primer ratio in seq. reaction



Comparison between Sequencing and Fragment Analysis Consumables



Application	Polymer	Capillary Length	Run Time	Throughput (24 hrs)		Resolution	Performance
				3100-Avant	3100		
High throughput, small size fragment analysis	POP-4™	22 cm	20 min	5,760* GT	23,040* GT	400 bp	0.50 SD†
Standard fragment analysis			45 min	2,560* GT	10,240* GT	400 bp	0.15 SD
Ultra rapid sequencing			40 min	72,000 bp	288,000 bp	500 bp	98.5% base calling accuracy
Rapid sequencing	POP-6™		60 min	48,000 bp	192,000 bp	500 bp	98.5% base calling accuracy
Standard sequencing	POP-4™	36 cm	100 min	34,500 bp	138,000 bp	600 bp	Q20 (KB™ Basecaller)
Long fragment analysis			65 min	1,760* GT	7,040* GT	500 bp	0.15 SD
Standard sequencing			2.5 hrs	23,400 bp	93,600 bp	650 bp	98.5% base calling accuracy
Long fragment analysis	POP-6™	50 cm	90 min	1,200* GT	4,800* GT	500 bp	0.15 SD
Long read sequencing	POP-4™	80 cm	3.5 hrs	22,800 bp	91,200 bp	950 bp	98.5% base calling accuracy

兩大主要應用：

DNA 定序
(Sequencing)

ATGC

基因體計畫
核心實驗室
(研究機構、大學、醫院)
代客定序公司
生物醫學、農業、小型研究

基因體計畫
核心實驗室
Gene Hunting
親子鑑定服務(檢調單位)
SSCP
AFLP等



Genotyping
(GeneScan)


- De Novo Sequencing
- MicroSeq 16 S 微生物菌種鑑定
- HIV Genotyping system
- Sequencing Based HLA typing
- Mitochondrial Sequencing
- Methylation Study



AB Applied Biosystems

16S rRNA Structure & Primer Locations

- Present in all bacteria
- Slow molecular clock
- Used for taxonomy
- Publish new species
- In public data bases




MicroSeq Libraries



- **16S Bacteria Library Eubacteria Branch**
 - ◆ Pharmaceutical / Sterility Control / Human pathogens
 - Over 1,100 entries, over 1,000 type strains.
 - Near complete Bacillus, Staphylococcus & Coryneforms (40 Corynebacteria)
 - Very strong in Gram negative nonfermenters.
 - Opportunistic pathogens, and 450-500 Aerobes.
 - Add Mycoplasmas, Ureaplasmas
- **28S Fungi Library**
 - ◆ 500 Yeasts.
 - ◆ 1,000 Filamentous Fungi.
- **Typing Kits**
 - ◆ Methicillin Resistant Staph. aureus (MRSA) Library/
Protein A gene used for typing

MicroSeq® Workflow


AB Applied Biosystems




Pharma Sterility Control



Genentech Unknowns D1796, D1921 & D1921B = *Serratia marcescens*



mitoSEQr™ Resequencing System



Nuclear DNA

- 2 copies/cell
- Unique to sample
- Bi-parental inheritance
- Recombination

Mitochondrial DNA

- 100 - >1000 copies/cell
- May not be unique to sample (heteroplasmy)
- Predominantly maternal inheritance
- No recombination
- High Mutation Rate: 10-20x nuclear DNA
- Highly polymorphic in Control Region

Hyper-variable regions
HV1 HV2

Control region

mtDNA

16,569 bases

Mitochondrial DNA Sequencing is used for

- Human Evolution & Origins
- Human mtDNA Disease Identification
- SNP and Mutation Detection
- Senescence
- Matching mtDNA sequence with maternal relatives – (even many generations apart)

Date | Confidential

© 2005 Applied Biosystems

Mitochondrial DNA SNP Detection

- One of the major applications is for cancer detection
 - ◆ mtDNA mutation occurs in large number and high frequency in cancers
 - ◆ Mutation not restricted by cancer type (*Hum Mutat. 2006 May;27(6):575-582*)
 - ◆ High copy number
 - ◆ Used for early cancer detection because mtDNA mutation occurs early in tumor progression (*J Mol Diagn. 2005 May;7(2):258-67. Mitochondrial DNA as a Cancer Biomarker*)
- Studies on other mitochondrial diseases

Date | Confidential

© 2005 Applied Biosystems

Mitochondrial Diseases

Applied Biosystems

Human mtDNA (16.5 kb)

<img alt="Diagram of the Human mitochondrial DNA (mtDNA) molecule, which is a circular molecule of approximately 16.5 kb. The diagram shows various genes and their locations: ND1, ND2, ND3, ND4L/ND4, ND5, ND6, Cyt b, and ATPase 6/8. It also shows tRNAs labeled L1, L2, L3, L4, L5, L6, L7, L8, L9, L10, L11, L12, L13, L14, L15, L16, L17, L18, L19, L20, L21, L22, L23, L24, L25, L26, L27, L28, L29, L30, L31, L32, L33, L34, L35, L36, L37, L38, L39, L40, L41, L42, L43, L44, L45, L46, L47, L48, L49, L50, L51, L52, L53, L54, L55, L56, L57, L58, L59, L60, L61, L62, L63, L64, L65, L66, L67, L68, L69, L70, L71, L72, L73, L74, L75, L76, L77, L78, L79, L80, L81, L82, L83, L84, L85, L86, L87, L88, L89, L90, L91, L92, L93, L94, L95, L96, L97, L98, L99, L100, L101, L102, L103, L104, L105, L106, L107, L108, L109, L110, L111, L112, L113, L114, L115, L116, L117, L118, L119, L120, L121, L122, L123, L124, L125, L126, L127, L128, L129, L130, L131, L132, L133, L134, L135, L136, L137, L138, L139, L140, L141, L142, L143, L144, L145, L146, L147, L148, L149, L150, L151, L152, L153, L154, L155, L156, L157, L158, L159, L160, L161, L162, L163, L164, L165, L166, L167, L168, L169, L170, L171, L172, L173, L174, L175, L176, L177, L178, L179, L180, L181, L182, L183, L184, L185, L186, L187, L188, L189, L190, L191, L192, L193, L194, L195, L196, L197, L198, L199, L200, L201, L202, L203, L204, L205, L206, L207, L208, L209, L210, L211, L212, L213, L214, L215, L216, L217, L218, L219, L220, L221, L222, L223, L224, L225, L226, L227, L228, L229, L230, L231, L232, L233, L234, L235, L236, L237, L238, L239, L240, L241, L242, L243, L244, L245, L246, L247, L248, L249, L250, L251, L252, L253, L254, L255, L256, L257, L258, L259, L250, L251, L252, L253, L254, L255, L256, L257, L258, L259, L260, L261, L262, L263, L264, L265, L266, L267, L268, L269, L2610, L2611, L2612, L2613, L2614, L2615, L2616, L2617, L2618, L2619, L2620, L2621, L2622, L2623, L2624, L2625, L2626, L2627, L2628, L2629, L2630, L2631, L2632, L2633, L2634, L2635, L2636, L2637, L2638, L2639, L2640, L2641, L2642, L2643, L2644, L2645, L2646, L2647, L2648, L2649, L2650, L2651, L2652, L2653, L2654, L2655, L2656, L2657, L2658, L2659, L2660, L2661, L2662, L2663, L2664, L2665, L2666, L2667, L2668, L2669, L2670, L2671, L2672, L2673, L2674, L2675, L2676, L2677, L2678, L2679, L2680, L2681, L2682, L2683, L2684, L2685, L2686, L2687, L2688, L2689, L2690, L2691, L2692, L2693, L2694, L2695, L2696, L2697, L2698, L2699, L2610, L2611, L2612, L2613, L2614, L2615, L2616, L2617, L2618, L2619, L2620, L2621, L2622, L2623, L2624, L2625, L2626, L2627, L2628, L2629, L26210, L26211, L26212, L26213, L26214, L26215, L26216, L26217, L26218, L26219, L26220, L26221, L26222, L26223, L26224, L26225, L26226, L26227, L26228, L26229, L26230, L26231, L26232, L26233, L26234, L26235, L26236, L26237, L26238, L26239, L26240, L26241, L26242, L26243, L26244, L26245, L26246, L26247, L26248, L26249, L26250, L26251, L26252, L26253, L26254, L26255, L26256, L26257, L26258, L26259, L26260, L26261, L26262, L26263, L26264, L26265, L26266, L26267, L26268, L26269, L26270, L26271, L26272, L26273, L26274, L26275, L26276, L26277, L26278, L26279, L26280, L26281, L26282, L26283, L26284, L26285, L26286, L26287, L26288, L26289, L26290, L26291, L26292, L26293, L26294, L26295, L26296, L26297, L26298, L26299, L26210, L26211, L26212, L26213, L26214, L26215, L26216, L26217, L26218, L26219, L26220, L26221, L26222, L26223, L26224, L26225, L26226, L26227, L26228, L26229, L26230, L26231, L26232, L26233, L26234, L26235, L26236, L26237, L26238, L26239, L26240, L26241, L26242, L26243, L26244, L26245, L26246, L26247, L26248, L26249, L26250, L26251, L26252, L26253, L26254, L26255, L26256, L26257, L26258, L26259, L26260, L26261, L26262, L26263, L26264, L26265, L26266, L26267, L26268, L26269, L26270, L26271, L26272, L26273, L26274, L26275, L26276, L26277, L26278, L26279, L26280, L26281, L26282, L26283, L26284, L26285, L26286, L26287, L26288, L26289, L262810, L262811, L262812, L262813, L262814, L262815, L262816, L262817, L262818, L262819, L262820, L262821, L262822, L262823, L262824, L262825, L262826, L262827, L262828, L262829, L2628210, L2628211, L2628212, L2628213, L2628214, L2628215, L2628216, L2628217, L2628218, L2628219, L2628220, L2628221, L2628222, L2628223, L2628224, L2628225, L2628226, L2628227, L2628228, L2628229, L26282210, L26282211, L26282212, L26282213, L26282214, L26282215, L26282216, L26282217, L26282218, L26282219, L26282220, L26282221, L26282222, L26282223, L26282224, L26282225, L26282226, L26282227, L26282228, L26282229, L262822210, L262822211, L262822212, L262822213, L262822214, L262822215, L262822216, L262822217, L262822218, L262822219, L262822220, L262822221, L262822222, L262822223, L262822224, L262822225, L262822226, L262822227, L262822228, L262822229, L2628222210, L2628222211, L2628222212, L2628222213, L2628222214, L2628222215, L2628222216, L2628222217, L2628222218, L2628222219, L2628222220, L2628222221, L2628222222, L2628222223, L2628222224, L2628222225, L2628222226, L2628222227, L2628222228, L2628222229, L26282222210, L26282222211, L26282222212, L26282222213, L26282222214, L26282222215, L26282222216, L26282222217, L26282222218, L26282222219, L26282222220, L26282222221, L26282222222, L26282222223, L26282222224, L26282222225, L26282222226, L26282222227, L26282222228, L26282222229, L262822222210, L262822222211, L262822222212, L262822222213, L262822222214, L262822222215, L262822222216, L262822222217, L262822222218, L262822222219, L262822222220, L262822222221, L262822222222, L262822222223, L262822222224, L262822222225, L262822222226, L262822222227, L262822222228, L262822222229, L2628222222210, L2628222222211, L2628222222212, L2628222222213, L2628222222214, L2628222222215, L2628222222216, L2628222222217, L2628222222218, L2628222222219, L2628222222220, L2628222222221, L2628222222222, L2628222222223, L2628222222224, L2628222222225, L2628222222226, L2628222222227, L2628222222228, L2628222222229, L26282222222210, L26282222222211, L26282222222212, L26282222222213, L26282222222214, L26282222222215, L26282222222216, L26282222222217, L26282222222218, L26282222222219, L26282222222220, L26282222222221, L26282222222222, L26282222222223, L26282222222224, L26282222222225, L26282222222226, L26282222222227, L26282222222228, L26282222222229, L262822222222210, L262822222222211, L262822222222212, L262822222222213, L262822222222214, L262822222222215, L262822222222216, L262822222222217, L262822222222218, L262822222222219, L262822222222220, L262822222222221, L262822222222222, L262822222222223, L262822222222224, L262822222222225, L262822222222226, L262822222222227, L262822222222228, L262822222222229, L2628222222222210, L2628222222222211, L2628222222222212, L2628222222222213, L2628222222222214, L2628222222222215, L2628222222222216, L2628222222222217, L2628222222222218, L2628222222222219, L2628222222222220, L2628222222222221, L2628222222222222, L2628222222222223, L2628222222222224, L2628222222222225, L2628222222222226, L2628222222222227, L2628222222222228, L2628222222222229, L26282222222222210, L26282222222222211, L26282222222222212, L26282222222222213, L26282222222222214, L26282222222222215, L26282222222222216, L26282222222222217, L26282222222222218, L26282222222222219, L26282222222222220, L26282222222222221, L26282222222222222, L26282222222222223, L26282222222222224, L26282222222222225, L26282222222222226, L26282222222222227, L26282222222222228, L26282222222222229, L262822222222222210, L262822222222222211, L262822222222222212, L262822222222222213, L262822222222222214, L262822222222222215, L262822222222222216, L262822222222222217, L262822222222222218, L262822222222222219, L262822222222222220, L262822222222222221, L262822222222222222, L262822222222222223, L262822222222222224, L262822222222222225, L262822222222222226, L262822222222222227, L262822222222222228, L262822222222222229, L2628222222222222210, L2628222222222222211, L2628222222222222212, L2628222222222222213, L2628222222222222214, L2628222222222222215, L2628222222222222216, L2628222222222222217, L2628222222222222218, L2628222222222222219, L2628222222222222220, L2628222222222222221, L2628222222222222222, L2628222222222222223, L2628222222222222224, L2628222222222222225, L2628222222222222226, L2628222222222222227, L2628222222222222228, L2628222222222222229, L26282222222222222210, L26282222222222222211, L26282222222222222212, L26282222222222222213, L26282222222222222214, L26282222222222222215, L26282222222222222216, L26282222222222222217, L26282222222222222218, L26282222222222222219, L26282222222222222220, L26282222222222222221, L26282222222222222222, L26282222222222222223, L26282222222222222224, L26282222222222222225, L26282222222222222226, L26282222222222222227, L26282222222222222228, L26282222222222222229, L262822222222222222210, L262822222222222222211, L262822222222222222212, L262822222222222222213, L262822222222222222214, L262822222222222222215, L262822222222222222216, L262822222222222222217, L262822222222222222218, L262822222222222222219, L262822222222222222220, L262822222222222222221, L262822222222222222222, L262822222222222222223, L262822222222222222224, L262822222222222222225, L262822222222222222226, L262822222222222222227, L262822222222222222228, L262822222222222222229, L2628222222222222222210, L2628222222222222222211, L2628222222222222222212, L2628222222222222222213, L2628222222222222222214, L2628222222222222222215, L2628222222222222222216, L2628222222222222222217, L2628222222222222222218, L2628222222222222222219, L2628222222222222222220, L2628222222222222222221, L2628222222222222222222, L2628222222222222222223, L2628222222222222222224, L2628222222222222222225, L2628222222222222222226, L2628222222222222222227, L2628222222222222222228, L2628222222222222222229, L26282222222222222222210, L26282222222222222222211, L26282222222222222222212, L26282222222222222222213, L26282222222222222222214, L26282222222222222222215, L26282222222222222222216, L26282222222222222222217, L26282222222222222222218, L26282222222222222222219, L26282222222222222222220, L26282222222222222222221, L26282222222222222222222, L26282222222222222222223, L26282222222222222222224, L26282222222222222222225, L26282222222222222222226, L26282222222222222222227, L26282222222222222222228, L26282222222222222222229, L262822222222222222222210, L262822222222222222222211, L262822222222222222222212, L262822222222222222222213, L262822222222222222222214, L262822222222222222222215, L262822222222222222222216, L262822222222222222222217, L262822222222222222222218, L262822222222222222222219, L262822222222222222222220, L262822222222222222222221, L262822222222222222222222, L262822222222222222222223, L262822222222222222222224, L262822222222222222222225, L262822222222222222222226, L262822222222222222222227, L262822222222222222222228, L262822222222222222222229, L2628222222222222222222210, L2628222222222222222222211, L2628222222222222222222212, L2628222222222222222222213, L2628222222222222222222214, L2628222222222222222222215, L2628222222222222222222216, L2628222222222222222222217, L2628222222222222222222218, L2628222222222222222222219, L2628222222222222222222220, L2628222222222222222222221, L2628222222222222222222222, L2628222222222222222222223, L2628222222222222222222224, L2628222222222222222222225, L2628222222222222222222226, L2628222222222222222222227, L2628222222222222222222228, L2628222222222222222222229, L26282222222222222222222210, L26282222222222222222222211, L26282222222222222222222212, L26282222222222222222222213, L26282222222222222222222214, L26282222222222222222222215, L26282222222222222222222216, L26282222222222222222222217, L26282222222222222222222218, L26282222222222222222222219, L26282222222222222222222220, L26282222222222222222222221, L26282222222222222222222222, L26282222222222222222222223, L26282222222222222222222224, L26282222222222222222222225, L26282222222222222222222226, L26282222222222222222222227, L26282222222222222222222228, L26282222222222222222222229, L262822222222222222222222210, L262822222222222222222222211, L262822222222222222222222212, L262822222222222222222222213, L262822222222222222222222214, L262822222222222222222222215, L262822222222222222222222216, L262822222222222222222222217, L262822222222222222222222218, L262822222222222222222222219, L262822222222222222222222220, L26282222222222222

Methylation Study



- Bisulfite treatment of gDNA converts unmethylated Cytosine to Uracil.
 - ◆ Methylated Cytosine remains a **C**.
- After PCR amplification all Uracils are converted Thymine.
- Comparing the sequences of the native gDNA to the bisulfite treated gDNA will show you the methylated bases.

	unmethylated	methylated
Target	-- <u>C</u> G--	-- <u>C^m</u> G--
Bisulfite treatment & cleanup	-- <u>U</u> G--	-- <u>C^m</u> G--
PCR	-- <u>T</u> G--	-- <u>C</u> G--


Date | Confidential

© 2005 Applied Biosystems

3100 Principles and Concepts

GENESCAN


Standards




Sample



3100 Principles and Concepts




45 min /each run



GeneMapper™ Analysis Software

AB Applied Biosystems


- Automatic size determination
- Automatic allele peak quantitation
- Used in conjunction with Genotyper™




Genotyping Applications 基因比對之應用

AB Applied Biosystems

➤ Human Linkage Mapping
➤ microsatellite genome panel



➤ 人類身分鑑定
➤ 16 STR Markers in one lane




➤ AFLP Plant Mapping Kits



➤ SNPs Detection

➤ SSCP

Genotype with Microsatellite Markers



Microsatellites are 2–7 bp repeats flanked with unique sequences constituting a 'microsatellite marker'


- » The number of repeat units varies in the population
- » Assumption is that each repeat unit is randomly distributed in the population

Widely spread in the genome

Over 7000 microsatellite markers have been mapped in the human genome



3130 GeneScan Analysis for Microsatellites



GeneScan Analysis plots (size (bp) vs. fluorescence intensity) of **microsatellites** on the 3100 Genetic Analyzer. GeneScan Installation Standard and the GeneScan 400-HD Size Standard run with POP4 and a 36cm capillary array are shown. The Installation Standard is a pool of six Linkage Mapping Set markers, amplified from an individual CEPH family member (1347-02).


Possible Mechanism for "Stuttering"




Step 1




Step 2



Step 3




Step 4



Date | Confidential


© 2005 Applied Biosystems

STR Stutter Peaks



Date | Confidential

© 2005 Applied Biosystems



AmpF[®]STR[®] Identifier[™] Kit A Combination of Experience

親子鑑定


SGM Plus [™] Kit	Profiler Plus [™] Kit	COfiler [™] Kit
D2S1338		
D16S539	FGA	D16S539
vWA	vWA	D3S1358
D3S1358	D3S1358	
D18S51	D18S51	CSF1PO
D21S11	D21S11	TPOX
D8S1179	D8S1179	TH01
		D7S820
FGA	D7S820	
TH01	D13S317	
D19S433	D5S818	



16 markers in
one PCR reactions




AB Applied Biosystems

AmpF[®]STR[®] Identifier[™] Allelic Ladder







Advantages of AFLP

Applied Biosystems

Some of the advantages of the AFLP technique are the following:

- ◆ Only small amounts of DNA are needed.
- ◆ Unlike randomly amplified polymorphic DNAs (RAPDs) that use multiple, arbitrary primers and lead to unreliable results, the AFLP technique uses only two primers and gives reproducible results.
- ◆ Many restriction fragment subsets can be amplified by changing the nucleotide extensions on the adaptor sequences. Hundreds of markers can be generated reliably.
- ◆ High resolution is obtained because of the stringent PCR conditions.
- ◆ The AFLP technique works on a variety of genomic DNA samples.
- ◆ No prior knowledge of the genomic sequence is required.

Template preparation & ligation of AFLP adaptor

A. Cut genomic DNA into fragments with the restriction enzymes Msel and EcoRI:




B. Ligate adaptors: EcoRI  and Msel 

C. Modify genomic DNA fragments:



Preselective amplification

Prepared Template: Genomic DNA Fragment, Modified with Adaptors



Preselective Primers:

 A : EcoRI adaptor + recognition site + A
or EcoRI adaptor + recognition site

 C : Msel adaptor + recognition site + C


Selective amplification

AB Applied Biosystems


A. Choose Selective AFLP Primers: ★ ■ Axx □ Cxx ★ Fluorescent dye

★ ■ Axx - one of sixteen different fluorescent dye-labeled AFLP EcoRI Selective Amplification primers.
□ Cxx - one of eight different AFLP MseI Selective Amplification primers.

B. Run Selective Amplification:



AB Applied Biosystems





AB Applied
Biosystems

Thanks!

www.appliedbiosystems.com.tw

twsupport@appliedbiosystems.com

0800-251326 ; 02-23582838