

ELISA的原理及衍伸應用及 Western Blot的常見問答集

洪國彥 產品專員

劉昌維 製程經理

亞諾法生技

什麼是ELISA?

- ▶ Enzyme-linked immunosorbent assay 酶素免疫分析法
- ▶ 利用固體承載物，以及抗體抗原專一性，對於樣品中的目標物進行檢測
- ▶ 利用最後呈色的量級進行定量



- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting WB 疑難排除


ELISA types (ELISA 之分類)

- ▶ Indirect (間接式) ELISA
- ▶ Direct (直接式) ELISA
- ▶ Sandwich (雙抗夾心) ELISA
- ▶ Competitive (競爭性) ELISA



Indirect (間接式) ELISA

- ▶ Antigen: immobilized on the plate
- ▶ Primary Antibody (一抗): IgG
- ▶ Secondary Antibody (二抗): Anti-IgG
- ▶ Indicator enzyme: HRP (Horse Radish Peroxidase)
- ▶ Chromogen: TMB




Sandwich (雙抗夾心) ELISA

▶ Antibody Pair

- Recognize different Epitopes (抗原表位): 不適用分子量較小抗原表位少之標的偵測物。
- Capture Antibody (用於捕獲抗原): 常用單株抗體
- Detection Antibody (用於偵測被捕獲之抗原): enzyme-labeled


▶ HRP (Indicator)

▶ TMB (Chromogen)



Competitive (競爭型) ELISA

- ▶ Capture Antibody/Antigen (捕獲)
- ▶ Competitor Required (競爭者):
 - ▶ antigen vs. Indicator-labeled antigen
- ▶ Color development: HRP & TMB



- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting WB 疑難排除

ELISA 內容物介紹

- ▶ Antibody/antigen coating (96-well) plate
- ▶ Standards/calibrators (標準樣本)
- ▶ Diluent (稀釋液)
- ▶ Capture & detection antibody
- ▶ Chromogen (顯色劑)
- ▶ Stop solution (終止液)




Sample Pretreatment (樣本前處理)

- ▶ 細胞培養液上清液, 血清
以離心去除非水溶性顆粒.
- ▶ 血漿
加入 citrate, EDTA 或 heparin 作為抗凝劑, 再離心並收集上清液.
- ▶ 檢測前需再經過稀釋
- ▶ 樣本經前處理後應盡速檢測, 或分裝後保存於-20°C.

Standard Preparation (標準樣本)

- 以稀釋緩衝液將標準樣本進行系列稀釋



ELISA常規步驟 (雙抗夾心)



樣本前處理

標準樣本
系列稀釋

包被

清洗

封閉



清洗

加入樣本/
標準品



清洗

加入檢測抗體

清洗


呈色

吸光度測定



Blocking (封閉)

- ▶ 封閉非專一性蛋白質結合位點
- ▶ Blocking Reagent:
 - 1% BSA (Immunoassay Grade)
 - 10% FBS
 - 5% skim milk in PBS
 - 5% serum



ELISA常規步驟 (雙抗夾心)



樣本前處理

標準樣本
系列稀釋

包被

清洗

封閉



清洗

加入樣本/
標準品



清洗

加入檢測抗體

清洗

呈色

吸光度測定



- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting WB 疑難排除



Data Analysis


- ▶ Standard curve (標準曲線)
- ▶ Precision (精確度)
- ▶ Accuracy (準確度)
- ▶ Sensitivity (靈敏度)
- ▶ Specificity (特異性)



Standard Curve (標準曲線)

	LipoDox conc. (ng/ml)							
Con.	0	3.13	6.25	12.5	25	50	100	
OD1	0.12	0.22	0.34	0.48	0.85	1.32	1.99	
OD2	0.12	0.19	0.29	0.55	0.83	1.17	1.84	
AVG	0.12	0.21	0.31	0.51	0.84	1.25	1.92	
CV	1%	10%	11%	9%	2%	9%	6%	

變異係數
(coefficient of variation)
Sample CV < 15%
Blank CV < 20%



Precision (精確度)

- ▶ Intra-assay (內部檢測) Precision (Precision within an assay)
- ▶ Inter-assay (批間檢測) Precision (Precision between assays)

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	4	5	6
n	20	20	20	60	60	60
Mean (nmol/L)	6.9	18.2	40.2	10.9	25.7	67.2
Standard deviation	0.3	1.3	2.5	0.7	1.4	3.8
CV (%)	4.3	7.1	6.2	6.4	5.4	5.7

Less than 10%

Accuracy 准確度(Recovery Studies)

- ▶ 樣本成分會影響吸亮度之測定
- ▶ Recovery: 80-120%

Sample	Average % Recovery	Range
Cell culture media	97	85 - 104%
Serum	98	93 - 103%
EDTA plasma	100	94 - 105%
Heaprin plasma	98	92 - 104%
Citrate plasma	101	92 - 106%



Sensitivity (靈敏度)

▶ Limit of Detection (檢測限)

樣品中所含的目標蛋白可被檢測出之**最低量或濃度**.



Specificity (專一性)

Metanephine (甲基腎上腺素) ELISA Kit (KA1890)

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Metanephine
	Derivatized Metanephine	100
	Derivatized Normetanephrine	0.15
	Derivatized 3-methoxytyramine	< 0.001
	Adrenaline	3.3
	Noradrenaline	< 0.001
	Dopamine	< 0.001
	Vanillic mandelic acid, L-Dopa, Homovanillic acid, L-Tyrosin, Tyramin	< 0.001



Abnova

www.abnova.com

- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting WB 疑難排除



ELISA 之應用

Application

- ▶ ELISPOT
- ▶ Cell-Based ELISA Kit

Related product in Abnova

- ▶ Ab Pair for ELISA (1000+)
- ▶ ELISA Kit (2000+)
- ▶ Cell-Based ELISA Kit (500+)




Abnova

www.abnova.com

ELISPOT (Enzyme-linked Immunospot Assay) (酶聯免疫斑點法)

- Legend**
- Y Antibody
 - Secreted Analyte
 - B Biotinylated Antibody
 - ◆ Alkaline phosphatase Conjugated Streptavidin
 - Color Product
 - BCIP/NBT



適合檢測外泌性蛋白


Cell-Based ELISA

將細胞接種到各孔的底部上。將細胞固定和封閉。

加入一抗和胞內之專一性抗原結合。

加入二抗。

HRP 酶將加入之 TMB 氧化。顏色由藍變黃，加入酸將反應終止，再測量吸光值 (450nm)。



Advantages of Cell-based ELISA Kit

	Cell-Based ELISA	Western Blot
樣本類型	活細胞即可 (Suspension cells, loosely attached cells, attached cells)	細胞裂解液
樣本前處理	不須額外處理	須將細胞裂解
檢測過程	快速	較費時 (Gel electrophoresis, Transfer)
檢測限	ng/mL	mg/mL
應用	Qualitative/phosphory- -lation detection	Qualitative/phosphoryl- ation detection



- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting WB 疑難排除



Weak or no signal 無信號

Possible Cause	Solution
The protein amount is added into well insufficiently.	Ensure extract contains enough amount of proteins.
Incubation time and temperature is incorrect.	Ensure the incubation time and temperature described in the protocol are correctly followed.
Inadequate reagent volumes or improper dilution	Check pipettes and ensure correct preparation



High background 高背景

Possible Cause	Solution
Poor blocking	Try a different blocking agent.
Contaminated wash buffer	Make fresh wash buffer
Insufficient washing	Increase number of washes Increase time of soaking between washes
Incubation time and temperature is incorrect.	Ensure the incubation time and temperature described in the protocol are correctly followed.
Primary and secondary antibodies	Try decreasing the antibody concentration and/or the incubation period.



Anomalous signal 讀值異常

Possible Cause	Solution
Reagents not fresh or not at the correct pH	Use freshly prepared solutions and reagents.
Reagent added in incorrect order, or incorrectly prepared	Review protocol and check the concentration of reagents.
Bubbles in the wells	Make sure you remove air bubbles within wells.
Uneven temperature around work surface	Avoid incubating plate in areas where environmental conditions vary, and use plate sealer
Filth at the bottom of 96-well plate	Make sure the surface is clean before adding the reagents into wells



- ELISA types 分類及原理
- ELISA Assay procedure 內容物與檢測介紹
- Data analysis 數據分析
- ELISA application ELISA 之應用
- ELISA trouble shooting 疑難排除
- Western Blot trouble shooting 疑難排除



什麼是Western Blot?

- 利用膠體電泳法分離樣品中的蛋白質、轉漬到膜上後，使用對目標蛋白具專一性的抗體進行蛋白質偵測。




Western Blot 起源?


- ▶ 史丹佛大學的George Stark團隊提出並發表於1979年。
 - They used passive transfer of the proteins, then ^{125}I -labeled Protein A for detection.
- ▶ 瑞士Harry Towbin團隊1979年發表接近於當代流程的方法
 - developed the method of electrophoretic transfer of proteins to membranes, their procedure also used secondary antibodies for detection.
- ▶ 西雅圖W. Neal Burnette
 - gave the technique the name “Western blotting” in 1979



Western Blot Procedure




Sodium Dodecyl Sulphate (SDS)



- SDS-PAGE 系統中，樣本分子的泳動率，
僅取決於其分子量，而與原來分子所帶的電荷無關

Sodium Dodecyl Sulphate (SDS)



SDS 分子上的非極性區 (綠色尾巴)，會與蛋白質上的非極性區結合，並且使蛋白質變性，成為一線狀長條分子，上面佈滿 SDS 分子。SDS 分子的極性區 (紫色圓點)，則露在外面，以增加親水性。


SDS-PAGE 系統中，樣本分子的泳動率，**僅取決於其分子量**，而與原來分子所帶的電荷無關。

Electrophoresis (電泳)

帶電分子在電場作用下，向著與其電性相反的電極移動，稱為電泳

$$\text{泳動率 } (\mathbf{v}) \sim \frac{\text{外加電壓 } (\mathbf{E}) \times \text{分子淨電荷 } (\mathbf{q})}{\text{分子泳動時與介質間之摩擦力 } (\mathbf{f})}$$

PAGE




利用一個分子篩來分離大小不同的蛋白質

Blocking (封閉)

- ▶ 封閉膜上非專一性蛋白質結合位點
- ▶ 減少抗體非專一性結合，進而減低背景干擾
- ▶ Blocking Reagent:
 - 1% BSA (Immunoassay Grade)
 - 10% FBS
 - 5% skim milk in PBS
 - 5% serum






1. Western blotting : 0.5ug /lane (12/18)						
Marker	1	2	3	4	5	6
Ag	GST	nv84	419	mouse serum	n98	nv8
1st Ab	Rat anti-GST(3000 X)					
2nd Ab	HRP-Goat-anti-Rat IgG(10000x)					

WB Troubleshooting

- ▶ Weak or no signal 無信號
- ▶ High background 高背景
- ▶ Multiple bands 多帶現象
- ▶ Speckled or Dotted Appearance 出現不均勻的斑點
- ▶ White bands on a black blot 出現白色條帶
- ▶ Smile effect of the bands 條帶 “微笑” 效應
- ▶ Bands Smeared, Bands not Sharp 條帶不銳利清晰



► Weak or No Signal 無信號





Troubleshooting

► Weak or No Signal 無信號

Possible Cause	Solution
Insufficient antigen	Load at least 20–30 µg protein per lane; Use protease inhibitors.
The 1 st Ab and the 2 nd Ab are not compatible	Use 2 nd Ab that was raised against the species in which the 1 st Ab was raised.
2nd Ab was inhibited by sodium azide	Do not use sodium azide together with HRP-conjugated antibodies.
Detection system	Check that the detection reagents are being stored correctly and used as recommended.



▶ High Background 高背景




Troubleshooting

▶ High Background 高背景

Possible Cause	Solution
Poor blocking	Try a different blocking agent (5% non-fat dry milk, 3% BSA, or normal serum).
Insufficient washing	Washing should be thorough at each step.
Choice of membrane	Nitrocellulose membrane is considered to give less background than PVDF.



► Multiple Bands 多帶現象



Abnova

www.abnova.com



Troubleshooting

▶ Multiple Bands 多帶現象

Possible Cause	Solution
Post-translation modification (acetylation, methylation, phosphorylation, glycosylation)	Examine the literature and use an agent to de-phosphorylate, de-glycosylate.
SDS caused nonspecific binding	Wash blots after transfer. Do not use SDS during immunoassay procedure.
Primary and secondary antibodies	Try decreasing the antibody concentration and/or the incubation period.



▶ Speckled or Dotted Appearance 不均勻的斑點



Troubleshooting


▶ Speckled or Dotted Appearance 不均匀的斑點

Possible Cause	Solution
Black Spots	Filter or change a new blocking reagent to avoid the antibodies are binding to blocking reagent.
White Spots	Make sure you remove air bubbles when preparing the gel for transfer.



▶ Other question

Smile effect of the bands



Bands Smeared



Abnova

www.abnova.com

Troubleshooting

Problem	Possible Cause	Solution
Smile effect of the bands	Migration is too fast or too hot	Slow down the migration or run the gel in the cold room or on ice.
Bands Smeared, Bands not Sharp	Bands smeared due to hot gel, high voltage	Decrease voltage, and prepare new running buffer.





Gene Name, ID, Catalog #, Alias

[Search](#) | [Advanced](#)


[Login](#) | [My Account](#) | [Interest product](#) | [Order Status](#) | [Cart](#)

Product Browser

- [DNA & RNA](#)
 - [siRNA/Chimera](#)
 - [Proteins/Peptides](#)
 - [MaxPab® Ab](#)
 - [Monoclonal Ab](#)
 - [Polyclonal Ab](#)
 - [Recombinant Ab](#)
 - [Conjugate Ab](#)
 - [Secondary/Tag Ab](#)
 - [Beads/Dyes](#)
 - [Lysates/Slides](#)
 - [Ab Pairs/Kits](#)
 - [Ab & Tissue Arrays](#)
 - [IHC/Pathology Ab](#)
 - [In Situ Hybridization](#)
 - [Circulating Rare Cell](#)
 - [Systems & Automations](#)
 - [GMP/ASR/IVD](#)
- Last updated: 2017/5/21

[Quick Order](#) ([Tutorial](#))


Chronic Lymphocytic Leukemia Microenvironment




Product & Service Showcase

All


T & B Cell Lineage Biomarkers




Broad Spectrum of Flow Cytometry Abs



Epithelial-Mesenchymal Transition



Prostate Cancer Antibodies



CAR T Cells



Mab in Oncology 300 selections

Promotion

SIGN UP TO WIN



- [Abnova Easter Egg Hunting Quest](#)
- [Scientific Image Contest](#)
- [New Products](#)
- [Reward Bonus](#)
- [More Promotions](#)

AbVideo™



Cell Culture (Attached Cells)

1st March to 31st August 2017 EST

INSPIRE

\$50 gift voucher +
Grand Prize
Echo Family

Thanks for your attention

