

HUMAN HEALTH | ENVIRONMENTAL HEALTH

In Vivo Imaging IVIS Lumina XRMS



產品應用專員 葉恭誌

0800-898-178 http://www.jnhtech.com.tw





- Why *in vivo* imaging
- Principles of optical in vivo imaging
- IVIS[®] Lumina XRMS Introduction
- Pre-clinical *in vivo* Imaging Agents

Application





Together, Bio-discovery and IVIS platforms provide the world leading portfolio for translational disease research, drug discovery and clinical development

... integrating all stages of research workflow

Current vs. BPI methodology



Current Methodology = 24 animals over four treatment points



Biophotonic imaging (BPI) Methodology = the same 6 animals over four treatment points



 Continued monitor Low variation

• 3R (Replacement, Reduction, Refinement)

In Vivo Imaging Landscape







Lumina XRMS



Quantitative 2D bioluminescence, fluorescence and X-ray imaging

Quantum GX



High Speed (8s), low Xray dose and High Resolution micro-CT

G8 PET-CT



Benchtop type High sensitivity 3D imaging Low cost of ownership

SOLARIS



Open air fluorescence imaging for Preclinical research in large animal models









Tailored To Therapeutic Applications





For the Bette



Non-invasive

- Does not require subject to be euthanized
- High throughput
 - Take 3 mouse at one time
- Multi-function
 - Bioluminescence

 fluorescence and X-Ray
- Easy to operation



Detecting MMP and cathepsin activity with MMPSense 680 and ProSense 750EX in orthotopic 4T1-luc2 tumors.



8



Principles of optical in vivo imaging

Tissue Absorption



Tissue Is Not Transparent – Light Absorbance Depends on Wavelength

Green laser reflection



Blue light transmission

Blue light reflection

Overview

Green laser transmission



Bioluminescence

Fluorescence





Tissue Is Not Transparent – Light Absorbance Depends on Wavelength

-HbO2

Renilla

Bacterial

-CB Green

Firefly

Luciferase

(Luc/2)

CB Red

(Lux)

Luciferase

Luciferase Spectrum

Optimal Imaging

Window

700

T 1.2

1

0.8

0.6

0.4

0.2

800

Normalized Intensity



Florescence Spectrum

Overview

1000000

100000

10000

1000

100

400

500

600

Wavelength (nm)

Hemoglobin μ_a (cm⁻¹)



The IVIS® views the diffuse image on the surface of the subject



- Light traveling through tissue scatters many times creating a "fuzzy" image at the surface of the animal.
- The IVIS[®] views the diffuse image on the surface of the subject.



How an X-Ray Image is Acquired





- X-rays will be attenuated in tissue differently resulting in an image on the scintillator
- The CCD views the scintillator resulting in a planar X-ray image
- X-ray and Optical images have different path lengths. To correct this geometrical difference, the X-ray image is registered to the optical image





Close Up of Scintillator in Position



IVIS® Lumina XRMS Introduction





CCD camera

Image chamber



anesthesia system



Acquisition computer & monitor





- High sensitivity CCD for bioluminescence or fluorescence imaging
- High throughput with 12.5 cm field of view
- > 26 filters, wavelength ranges from 500 865 nm
- Spectral unmixing using discrete bandpass filters
- Ideal for imaging multiple probes/reporters
- Software controlled settings and analysis tools





Software

Image Wizard: step by step for imaging setup

🔽 MIS Acquisition Control Pan	el		Simaging Wizard	
Imaging Mode Exposure Time ✓ Luminescent ✓ Fluorescent ✓ Photograph ▲uto ♀ Standard-On	Binning F/Stop Medium Binning F/Stop Medium Binning F/Stop Medium Binning F/Stop Medium Binning F/Stop Binning F/Stop	Excitation Filter Block	<u>\$</u>	Draging Mode Bioluminescence Imaging Select this option for imaging bioluminescent or chemiuminescent reporters, such as frefly kulferase, click beetle kulferase, renila, or bacterial kulferase.
Voverlay Lights Field of View: C	Alignment Grid System Status Idle	X-Rays will be produced when energized. Acquire	Puorescence	Select this option for imaging fluorescent proteins, dyes, or nanoparticles in the wavelength range of 450-050 nm. Both opidiumnuton (dumanationfrom above) and
Service 13.4 cm Subject height: 1.50 🖨 cm Focus: use subject height 🔻	Temperature:	Imaging Wizard	Cherenkov	Cherenkov Imaging Select this option for imaging radiotracers which

Living Image 4.5.2 download:

http://www.perkinelmer.com/lab-products-and-services/resources/software-downloads.html



Photographic /X-Ray + Optical Image = Overlay









Controls Sensitivity

🜠 IVIS Acquisition Control Panel								
Imaging Mode	Exposure Time	Binning	F/St	op	Excitation Filter	· Emi	ission Filter	
📃 🗹 Luminescent	1.00 💲 sec	Medium	∨ 1	¥	Block	🔽 Op	en	*
Eluorescent								
🚺 🗹 Photograph	Auto 🤤	Medium	~ 8	*				
Structure								
🗹 Overlay 📃 Lights 🗹 Alignment Grid								
Field of View: C	~	5	stem !	Statu	5		Acquire	
Service 12.9 cm		Idle			Acquiro			
							Imaging Wi	zard
Subject height: 1.50	ı 🔷 cm						Sequence	Setup
Focus: use subject	: height 🔽	Temperature:			Locked		Initializ	8



Settings



Camera and Lens Settings are Analogous to Those Used in Standard Photography

- Field of View (FOV) is dependent on the distance from the lens to the sample
- Light collected is proportional to how long the shutter is open (exposure time)
- Aperture (*f*/stop) controls the amount of light collected
- Digital pixel binning is possible on the CCD – alters sensitivity/resolution



Field of View (X-Ray)







- ▶ The IVIS[®] CCD camera has a raw signal range of 0 to 65,535 Analog to Digital counts (2¹⁶)
- Adjust camera settings to obtain a signal level of <u>600 to 60,000 counts</u> to be within the linear range of the detector
- Settings that control signal level are:
 - Exposure time
 - Pixel binning (CCD resolution)
 - f/stop (aperture)
- Instrument is calibrated to automatically compensate for changes in sensitivity settings when count levels are within the linear range



Exposure Time

- Signal level is directly proportional to exposure time (1:1)
- Shorter exposure time improves throughput
- Recommended minimum exposure time > 0.5 seconds
- Longer exposure times increase signal intensity
- Recommended maximum exposure time < 5 minutes



🚺 IVIS Acquisition Control Panel . Imaging Mode Exposure Time Emission Filter Binning F/Stop Excitation Filter 🗌 🗹 Luminescent 1.00 🗘 sec 🗸 ~ Medium Block Open ¥ Eluorescent 🚹 🗹 Photograph 🛛 Auto 😂 Medium 🔽 8 ¥ Structure V Overlay 📃 Lights 🔽 Alignment Grid Field of View: C System Status Acquire Idle 12.9 Service cm. Imaging Wizard 😂 cm Subject height: 1.50 Sequence Setup

Locked

Initialize



Temperature:

Focus: use subject height 🗸

Exposure time setting



Pixel Binning (CCD Resolution)



Binning refers to the grouping of pixels into a larger super-pixel

Pixel binning setting

Large Binning (16) Higher Sensitivity/Lower Resolution Medium Binning (8) Small Binning (4) Lower Sensitivity/Higher Resolution













f/stop (Lens Aperture)



f/stop (aperture) setting

f/stop controls the amount of light received by the CCD detector

- f/1 is wide open, maximum light collection – default for luminescent
- f/8 is smallest aperture, best resolution default for photo

🜠 IVIS Acquisition Control P	Panel			_ 🛛			
Imaging Mode Exposure Time	Binning F/Stop	Exc	itation Filter	Emission Filter			
📃 🗹 Luminescent 🛛 1,00 🤤 sec	🔽 Medium 🔽 1 🔍	Block	< 🗸 🗸 🗸	Open 🔽			
Eluorescent	1 2						
🚺 🗹 Photograph 🛛 Auto 🤤	Medium 🔽 4						
Structure ✓ Overlay Lights ✓ Alignment Grid							
Field of View: C	System Status		Acquire				
MIS	Idle						
Service 12.8 cm			Imagi	ng Wizard			
Subject height: 1,50 🔷 cm			🛛 🏼 🗰 Sequ	ience Setup			
Focus: use subject height 💌	Temperature:	Locked	Ir	itialize			



f/8





26



Radiance = Emission Light (photons/sec/cm²/str)

- Living Image[®] automatically compensates for <u>device settings</u>: Exposure time, *f*/stop, binning and field of View.
- Calibrated units are Photons per Second, representing the flux radiating omni-directionally from a user-defined region.
- Absolute calibration calibrated to the NIST standards







Raw Signal (Counts)











Exp time: 30 sec 30 sec 60 sec 60 sec 60 sec 60 sec **Binning:** small small small small medium medium Day: 2 3 4 5 6 1

Peak **Counts**





Calibrated Signal (Photons per second)













Exp time: Binning: Day:

30 sec small 1 30 sec small 2 60 sec small

3

60 sec 60 sec small medium 4 5

60 sec medium 6

Radiance: Photons per second







Excitation and Emission Filters



>26 filters, wavelength ranges from 500 - 865 nm

Lumina III Instruments					
Excitation Filters		Emission Filters			
Center λ	Band pass	Center λ	Band pass		
420	20	520	40		
440	20	570	40		
460	20	620	40		
480	20	670	40		
500	20	710	40		
520	20	790	40		
540	20	845	40		
560	20				
580	20				
600	20				
620	20				
640	20				
660	20				
680	20				
700	20				
720	20				
740	20				
760	20				
780	20				



• 27 commonly used fluorescent probes









GFP Well Plate Uncorrected



VS.

Units of 'Radiant Efficiency' compensates for nonuniform excitation light pattern

GFP Well Plate Corrected







Fluorescence

Animal Diet Autofluorescence in Control Mice





Ex 640 nm / Em 680 nm Ex 745 nm / Em 800 nm



- 13 - 19 30⁹ - 19 30⁹ - 03 - 03 - 03

(p/sec/cm²/lat)

Fluorescence

Mouse Hair Can Interfere With Both BLI and FLI






Sensitivity is a Function of a Signal to Noise



Luminescent Sources: Signal brightness generally lower than fluorescent sources

Higher sensitivity due to low level noise: both instrument and animal autoluminescence

Fluorescence Sources: Signals generally brighter than luminescent sources

Lower sensitivity due to higher noise: instrument background and autofluorescence





300

pixels

400

500

nha

200

100

700 · 600 · 500 · 400 · 300 · 200 · 100 ·

Improvements to Signal to Noise Ratio



Adaptive FL Background Subtraction: Software tool to reduce instrument background Spectral Unmixing: Extracts fluorescent signal from autofluorescence









What is Spectral Unmixing?







Autofluorescence







XenoLight 680/750 in a Mouse Using spectral unmixing



Raw Spectral Images





 Subcutaneous injections of 10¹⁴ molecules of XenoLight 680 (scruff)

 Subcutaneous injection of 10¹⁴ molecules of XenoLight 750 (lower dorsal region)

▶605nm excitation filter



Mechanics

- 1. Choose reporters that maximize signal-to-noise (S:N) ratio
- 2. Consider the appropriate control groups and imaging time points necessary
- 3. Use hairless mice or white-furred animals and depilate or shave
- 4. Switch to autofluorescence-free mouse diet
- 5. Closely map the kinetics of your biological bioluminescent
- 6. Animal handling can significantly affect kinetics
- 7. Image in the animal orientation that yields the highest sign
- 8. Cover intense signal to allow dimmer signals to dictate auto-exposure
- 9. Utilize guards to prevent reflection off neighboring animals
- 10. Use calibration unit for quantitation (Radiance and Radiant Efficiency)







Overview



Pre-clinical in vivo Imaging Agents

Imaging Basics – Reporter Molecules











XenoLight D-Luciferin (buy one get one free!)

Product	Product Description	Catalog Number
XenoLight D-Luciferin (K+ Salt) (1-50 g)	Lypholized bioluminescence substrate for in vivo imaging with Firefly Luciferase, in bulk	122799

- The ONLY Luciferin optimized for *in vivo* image with IVIS system
- Used by PerkinElmer scientists to calibrated IVIS imaging system
- Validated on our Bioware and Bioware Ultra Cell lies
- Supplied under stringent lot control to ensure experimental reproducibility

RediFect lentiviral particles

Product	Product Description	Catalog Number
RediFect Red-Fluc-Puromycin	Lentivirus particles containing red-shifted firefly luciferase with puromycin as selection marker	CLS960002
RediFect Red-Fluc-GFP	Lentivirus particles containing red-shifted firefly luciferase and Green Fluorescent Protein (GFP)	CLS960003
RediFect Green Renilla-Puromycin	Lentivirus particles containing Green Renilla luciferase with puromycin as selection marker	CLS960004

Bioware - Oncology Cell Lines



Bioware® Brite cell lines labeled with enhanced Red-Fluc vector

Product	Product Description	Catalog Number
HT1080-Red-Fluc	Human Fibrosarcoma Cancer Cell line.	BW 128092
4T1-Red-Fluc	Murine Breast Cancer Cell line	BW 124087
GL261-Red-Fluc	Murine Glioma Cell line	BW 134246
HepG2-Red-Fluc	Human Hepatic Cancer cell line	BW 134280
PC-3-Red-Fluc	Human Prostate Cancer Cell line	BW 128444
LnCaP-Red-Fluc	Human Prostate Cancer Cell line	BW 125055
B16-F10-Red-Fluc	Murine Melanoma Cancer Cell line	BW 124734
HCT-116-Red-Fluc	Human Colorectal Cancer Cell line	BW 124318
HT-29-Red-Fluc	Human Colorectal Cancer Cell line	BW 124353
Colo205-Red-Fluc	Human Colorectal Cancer Cell line	BW 124317
U-87 MG-Red-Fluc	Human Brain Cancer Cell line, ideal for glioblastoma models	BW 124577
NCI-H460-Red-Fluc	Human Lung Cancer Cell line, ideal for orthotopic lung tumor models	BW 124316
K-562-Red-Fluc	Human Leukemia Cell line	BW 124735
BxPC3-Red-Fluc	Human Pancreatic Cancer Cell	Lanacous
MCF-7-IRed-Fluc	Human Breast Cancer	
A549-Red-Fluc	Human Lung Cancer	
LL/2-Red-Fluc	Murine Lung Cancer	- 10 e d0 ⁷
SKOV3-Red-Fluc	Human Ovarian Cancer	0

4T1-luc2 - Tumor metastasis post intracardiac injection of cells (IVIS Lumina XR), MicroCT image (Quanutm FX) confirming bone degradation in the right tibia.

Color Scale Mm = 0.00e7 Mm = 1.40e6

Bacteria



Bacterium	Parental strain	Catalog No.
E. coli	EPEC WS2572 (Xen14) ETEC WS2583 (Xen16)	119223 119225
L. monocytogenes	ATCC 23074 (Xen19) 10403S (Serotype 1/2a wild-type strain) (Xen32)	119237 119238
P. aeruginosa	ATCC 19660 (Xen5) PAO1 (Xen41)	119228 119229
P. mirabilis	ATCC 51286 (Xen44)	119236
S. aureus	8325-4 (Xen8.1) ATCC 12600 (Xen29) ATCC 33591 (Xen31) ATCC 49525 (Xen36) UAMS-1 (Xen40)	119239 119240 119242 119243 119244
S. dysenteriae	88A6205. Clinical isolate (Xen27)	119231

Bacterium	Parental strain	Catalog No.
S. pyogenes	Strain 591, Group A, Serotype M49 (Xen20)	119250
S. typhimurium	FDA1189 (Xen33)	119235
Y. enterocolitica	91A1854 Clinical isolate (Xen24) WS2589 (Xen25)	119232 119233



Xen05: Pseudomonas aeruginosa





Xen44: Monitoring migration Xen: of UTI infection from infec the bladder to the kidney non-invasively in real time

Xen5: P. aeruginosa infection on a biofilm



Agent Categories



Agent Platforms: Vascular Fluorescent Agents



9_ 9_	Agent	Agent Description
Accession Accession	Genhance™ 680/750	Genhance™ 680/750: Small molecule fluorescence agent. Use as a control or in vascular imaging
	Superhance™ 680	Superhance™ 680: Small molecule agent. Binds to albumin in blood for extended (30m–1h) circulation
Linker	AngioSense [®] 680EX/750EX	AngioSense [®] 680/750: Agent that remains localized in vasculature for 0–4h.; 5 days in arthritic joints
AngioSense®	AngioSPARK [®] 680/750	AngioSPARK [®] 680/750:Pegylated fluorescent nanoparticles; Remains localized in vasculature
	GastroSense™ 750	Agent to monitor gastric emptying and the impact of various drugs on gastric motility

Targeted Fluorescent Agents



Agent	Binds to …	
BombesinRSense 680	Bombesin receptors	•
HER2Sense 645	HER2/Neu receptor	
FolateRSense 680	Folate Receptor Protein	
TlectinSense 680	Vascular Endothelial cells (N-actelyglucosamines)	Ц
OsteoSense [®] 680/750/800	Hydroxyapatite	
IntegriSense™ 655/680/750	Integrin αvβ3 antagonist	
BacteriSense 645	Negatively charged phospholipids in Bacterial membrane	
Annexin-Vivo 750	Phosphatidylserine during early apoptosis	
HypoxiSense 680	Carbonic Anhydrase IX in hypoxic tissue and cells	
COX-2 Probe	Cyclooxygenase-2 (COX-2)	
2-DG 750	Glucose uptake Imaging	
Transferrin-vivo	Transferrin receptors	HE





HER2/Neu+ tumor targeting by HER2Sense 645



Agent	Agent Description
MMPSense™ 680	MMPSense [™] 680: Activated by matrix metalloproteinases including MMP's MMP-2, -3, -9 and -13
MMPSense™ 645, 750 FAST	MMPSense™645m 750 FAST (Fluorescent Activatable Sensor Technology) is an MMP activatable agent
ProSense [®] 680/750	ProSense [®] 680/750: Activated by proteases: cathespins B, L, S, and plasmin
Neutrophil Elastase 680 FAST™	FAST agent activated by elastase produced by neutrophil cells
Cat B 680/750 FAST™	Cathepsin B selective FAST activatable agent
Cat K 680 FAST™	Cat K 680 FAST (Fluorescent Activatable Sensor Technology) is a Cathepsin K activatable agent
ReninSense680 FAST™	ReninSense680 FAST™ (Fluorescent Activatable Sensor Technology) is a renin activatable agent
PSA 750 Fast	Activatable agent that detects active PSA in vivo



Monitor protease activity associated with disease state

In Vivo Agent Applications



	Oncology	
	Her2Sense	
	HypoxiSense	
Cardiovascular	IntegriSense	
	Annexin Vivo	
	FolateR-Sense	
	ProSense	Inflammation
	ProSense FAST	
	CatB FAST	
	MMPSense	
	MMPSense FAST	
	AngioSPARK	
	AngioSense	
	Neutrophil Elastase FAST	
	Cat K FAST	
Arthritis	OsteoSense	Bone Biology
Hypertension	ReninSense	
	GastroSense	Gastric emptying
Bacterial Infection	BacteriSense	



PerkinElmer offers four categories of fluorescent *IN VIVO* imaging agents:



VivoTag[™] 680XL Protein Labeling Kit : designed for preparing fluorescently labeled antibodies, proteins or peptides for small animal in vivo imaging applications.



VivoTrack 680 : cell labeling agent that intercalates into the plasma membrane of primary cells and cell lines.



https://www.addgene.org/fluorescent-proteins/in-vivo/

The nonprofit plasmid repository		Login Create Account		Search for plasmids Q		
F	ind Plasmids +	Deposit Plasmids +	How to Order -	Plasmid Reference	- About Add	lgene 🗸
Specia	Special Collections / Fluorescent Protein Guide / In Vivo Imaging Guide					
Fluorescent Protein Guide: In Vivo Imaging						

Near-Infrared

Protein	Excitation/Emission	Brightness	Find Plasmids
iRFP670	643/670	12.7	piRFP670-N1
iRFP682	663/682	10.2	piRFP682-N1
iRFP702	673/702	7.6	piRFP702-N1
iRFP (aka iRFP713)	690/713	6.2	piRFP
iRFP720	702/720	5.8	piRFP720-N1
iSplit	690/713	5.3	pPAS-E and pK-GAFm





HUMAN HEALTH | ENVIRONMENTAL HEALTH



In Vivo Imaging Applications

August 2, 2017

© 2009 PerkinElmer



Longitudinal Monitoring of Tumor Development and Metastasis

Longitudinal monitoring of tumor development



Bioware Ultra: 4T1-luc2



5 cells



With Bioware Ultra one can start collecting data from Day 0, while with caliper measurements one has to wait at least 28 days to see any tumor growth

Longitudinal Monitoring of Brain Tumors



MicroCT co-registration of BLI orthotopic Brain tumor signal by U87-MG-luc2 cells



- Monitor early tumor behavior non-invasively
- Brain tumor growth and therapeutic effect *in vivo*
- 3D quantitative tomography for accurate tumor size and location



Non-Invasive Detection of Micrometastasis





 5×10^5 4T1-luc2-1A4 cells orthotopically into the abdominal mammary fat pads

Imaging Breast Cancer Metastases





45 days following IC delivery of MDA-MB-231 cells into immune-deficient mice





Fluorescence Imaging of Drug Efficacy in a Pancreatic Tumor Model





Therapeutic effect of DHA (Dihydroartemisinin), a semi-synthetic derivative and active metabolite of artemisinin (a naturally occurring compound) on **RFP-labeled** pancreatic tumor cells (BxPc3-RFP) studies *in vitro* and *in vivo*



DHA inhibits cell and tumor growth by interfering with cell proliferation and inducing apoptosis.



Treatment Effects Assessment by optical imaging



- Use Angiosense to evaluate CT-322 Treatment
- Annexin-Vivo detects Therapy-Induced Tumor Apoptosis





Annexin-Vivo detects and quantifies tumor death in living mice after as little as a single dose of CY, before any tumor size changes occur



In Vivo Tumor Vasculature Targeted PET/NIRF Imaging



⁶⁴Cu-MSN-800CW-TRC105(Fab) Mesoporous silica nanoparticle

Scheme 1. Schematic Illustration of the Synthesis of ⁶⁴Cu-MSN-800CW-TRC105(Fab)^a





4T1 Murine Breast Cancer Model

Feng Chen et al Mol. Pharmaceutics 2014, 11, 4007 - 4014



Stem cell application











Inflammation

Degree of Roasting is the Main Determinant of the Effects of Coffee on NF-KB



Transgenic NF-kB reporter mice were given a single dose of coffee extract by oral gavage 3h prior to s.c. LPS injection (at 0h). A) The luciferase activity was measured by in vivo imaging at 0h, 3h and 6h.

Paur et al, 2010



Multiplex fluorescence Imaging



Vascular Leak



Dual imaging of infection (BLI) and immune response (FLI)





Non-invasive assessment of bioluminescent *S.aureas* infection and EGFP PMN infiltration

Dual Reporter: Bacterial luc and GFAP Brain Imaging



luc and *lux* in Pneumococcal Meningitis





Simultaneous visualization of:

- Streptococcus pneumonia (lux)
- Glial Fibrillary Acidic Protein FAP (luc)


Others

Imaging of stem cells and regeneration in adult zebrafish











Enterotoxigenic *E. coli* (ETEC)



- 1. Choose reporters that maximize signal-to-noise (S:N) ratio
- 2. Consider the appropriate control groups and imaging time points necessary
- 3. Use hairless mice or white-furred animals and depilate or shave
- 4. Switch to autofluorescence-free mouse diet
- 5. Closely map the kinetics of your biological bioluminescent
- 6. Animal handling can significantly affect kinetics
- 7. Image in the animal orientation that yields the highest sign
- 8. Cover intense signal to allow dimmer signals to dictate auto-exposure
- 9. Utilize guards to prevent reflection off neighboring animals
- 10. Use calibration unit for quantitation (Radiance and Radiant Efficiency)







Thank you for your attention!

